Instructions For Use

REF: LPA 001

AneuCyte™ Prenatal Enumeration Probe kit

FOR PROFESSIONAL USE ONLY


Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Specification

X centromere Xp11.1 - q11.1 (DXZ1) Green
Y centromere Yp11.1 ó q11.1 (DYZ3) Orange
18 centromere 18p11.1 ó q11.1 (D18Z1) Blue

13 unique sequence (13q14.2) Green
21 unique sequence (21q22.13) Orange

AneuCyte™ is a prenatal assay for the rapid detection of trisomy 13, 18, 21 present in Patau, Edward and Down syndrome respectively as well as sex chromosome aneuploidies (such as Klinefelter and Turner syndromes). The assay needs to be used in conjunction with foetal karyotype analysis. AneuCyte is designed for fluorescence *in situ* hybridisation of interphase nuclei of uncultured amniotic fluid cells and will provide a result within 24 hours of receiving an amniocentesis sample. This test does not detect structural chromosome abnormalities, mosaicism and numerical abnormalities of other chromosomes. Reporting and interpretation of FISH should be consistent with professional standards of practice. The test is not designed to

evaluate the risk of trisomy and therapeutic action should not be initiated on the basis of FISH result alone.

AneuCyte contains two probe sets for two separate hybridisations. The probe set X, Y and 18 DNA is a mixture of green, orange and blue directly labelled fluorescent DNA probes specific for the alpha satellite DNA sequences at the DXZ1, DYZ3 and D18Z1 regions of the chromosomes X, Y and 18 respectively. The probe set 13 and 21 is a mixture of green and orange directly labelled fluorescent DNA probes containing unique sequences. The chromosome 13 probe includes RB1 gene.

The probe sets are intended for the detection and quantification of chromosomes 13, 21, 18, X and Y by fluorescence *in situ* hybridisation.

Materials provided

Probe: 100µl per vial (10 tests) or 300µl per vial (30 tests)

The probe sets are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Probe Set 1 (Tube labelled with yellow dot):

DXZ1 (X centromere α-satellite) labelled with a Green fluorophore (FITC spectrum): 12-15 ng/test
DYZ3 (Y centromere α-satellite) labelled with an Orange fluorophore (TRITC/Orange spectrum): 8-10 ng/test
D18Z1 (18 centromere α-satellite) labelled with a Blue fluorophore (Aqua spectrum): 112-140 ng/test

Probe Set 2:

13q14.3 unique sequence (including the markers D13S915, D13S1155, D13S1195) labelled with a Green fluorophore (FITC spectrum): 120-150 ng/test
21q22.2 unique sequence (including the markers D21S270, D21S1867, D21S337, D21S1425, D21S1444, D21S341) labelled with an Orange fluorophore (TRITC/Orange spectrum): 40-50 ng/test

DAPI/Antifade ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole) in vectashield, 500 µl per vial

The Aquarius AneuCyte kit contains 2 vials of probe and 1 vial of counterstain.

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl ó 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendations

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The triple bandpass filter DAPI/FITC/TRITC is optimal for viewing orange and green fluorophores simultaneously. The blue fluorophore has specificity to the Aqua or DEAC spectrum (Aqua filter is required). The triple bandpass filter DAPI/FITC/Texas Red can be used to view all three fluorophores simultaneously.

Sample Suitability and collection

The AneuCyte X/Y/18 and 13/21 kit is designed for use on uncultured amniotic fluid samples fixed in Carnoy's fixative (see procedure below). Amniotic fluid sample collection should be performed according to the institution guidelines.

Amniotic fluid samples that appear bloody or brown should not be used, since they may contain maternal blood and may lead to false results.

Suggested Protocol

Preparation of fresh amniotic fluid samples for FISH:

- 1) Centrifuge 2-5 ml of whole amniotic fluid (AF) specimen for 7 minutes at 180xg, carefully remove the supernatant without disturbing the cell pellet.
- 2) Resuspend the pellet in 2ml of 0.075M KCl. Leave at room temperature for 5 minutes.
- 3) Add 2 ml of fresh fixative (3:1 methanol:glacial acetic acid) to the cells/hypotonic solution, adding the first ml dropwise whilst continuously mixing. Mix well.
- 4) Centrifuge the suspension for 5 minutes at 280xg, carefully remove the supernatant and resuspend the pellet in 2 ml fresh fixative.
- 5) Fixed specimens can be stored at this stage in a freezer at -20°C.
- 6) If the sample is not to be frozen, centrifuge tube at 280xg for 5 minutes. Remove as much supernatant as possible without disturbing the cell pellet. Flick the tube to resuspend the pellet in the small amount of fluid remaining.
- 7) To prepare slides for FISH, drop the cell suspension directly onto slide, making 2 hybridization areas. Allow to air dry.

For late gestational age specimens additional slide pretreatment may be beneficial:

- 1) Place slide(s) prepared from uncultured amniocytes, in 2x SSC for 1 hour at 37°C.
- 2) Place slide(s) in freshly made pepsin working solution (5 mg pepsin added to 100 ml of 0.01 M HCl) for 13 minutes at 37°C.
- 3) Rinse slide(s) in phosphate buffered saline (PBS) at room temperature for 5 minutes.
- 4) Place slide(s) in post fixation solution (0.95% formaldehyde: 1.0 ml of 37% formaldehyde, 0.18g MgCl₂ and 39.0 ml of PBS. Store at 4°C. Use within 1 month.) for 5 minutes at room temperature.
- 5) Rinse slide(s) in PBS at room temperature for 5 minutes.
- 6) Immerse slide(s) in 70% ethanol at room temperature. Allow the slide(s) to stand in the ethanol wash for 1 minute.
- 7) Remove the slide(s) from 70% ethanol. Repeat step 6 with 85% ethanol, followed by 100% ethanol.

FISH Protocol

Slide preparation:

- 1) Wash slide in 2xSSC for 2mins at room temperature.
- 2) Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation:

- Remove probe from the -20°C freezer and allow to warm up to RT.
- Ensure probe solution is uniform by repeated, gentle pipette mixing.
- Remove 10 µl of probe per test, place in a microcentrifuge tube and return the remaining probe to -20°C freezer.
- Place probe and sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 mins.
- Spot 10 µl of the probe mixture onto the cell sample and carefully apply 24 x 24 mm coverslip. Seal with rubber solution glue and allow to dry completely.

Denaturation:

- Denature sample and probe simultaneously by heating slide on a 75°C (+/-1°C) hotplate for 2 mins.

Hybridisation:

- Hybridise slide overnight in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes:

- Remove coverslip and all traces of glue carefully
- Wash slide in 0.4xSSC (pH 7.0) at 72°C (+/-1°C) for 2 mins.
- Drain slide and wash in 2xSSC, 0.05% Tween-20 (pH 7.0) at RT for 30 seconds
- Drain the slide and apply 10 µl of DAPI antifade
- Cover with a coverslip and allow colour to develop in the dark for 10 mins
- View with fluorescence microscope

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural recommendations

- Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Ltd.
- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Interpretation of Results

The sensitivity and specificity of FISH depends on a number of parameters, which vary from one cell type to another; from one probe to another; with the cell techniques used; and within the individual laboratory. Therefore we recommend for the use of AneuCyte that each laboratory must have its own standard material and determine its own FISH assay cutoff values for karyotypically normal and aneuploid samples (for guidelines contact Cytocell).

Customer Support

Please contact the Cytocell Sales and Marketing Department.

Patents and Trademarks

Aquarius and Cytocell are registered trademarks of Cytocell Ltd.

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