Amnio dish

Protocollo di utilizzo
Protocole d'utilisation
Benutzungsprotokoll
Protocol of use
Protocolo de utilización
Protocolo de utilização
Πρωτόκολλο χρήσης

- PER COLTURADI CELLULE
- POUR LA CULTURE DES CELLULES
- FÜR ZELLKULTUREN
- FOR CELL CULTURE
- PARA CULTIVO DE CÉLULAS
- PARA CULTURA DE CÉLULAS
- ΓΙΑ ΚΑΛΛΙΕΡΓΕΙΑ ΚΥΤΤΑΡΩΝ

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Sterile Petri dish with glass slide for culture and analysis of amniocytes and chorionic villi

Ref. EK AMN 240 240 tests

Material supplied

Each box contains 40 packages ready for use, for a total of 240 Amniodish.

Slide dimensions: 32 mm diameter, 0.13-0.16 mm thickness.

Validity: 2 years.

Amniocyte culture in situ protocol

Note: thaw the Amniomed Plus (ref. EKAMG200) in a 37°C water bath and mix well by swirling prior to use. Addition of L-Glutamine, antibiotics and serum are not necessary since these components are already present.

1. Place 0,5 ml of cell suspension gently onto each Amniodish, on the glass slide.
2. Incubate at 37°C, 5% CO₂ for 6 - 12 hours, with minimal mechanical disturbance.
3. Add 2 ml of fresh Amniomed Plus.
4. On Day 6, replace the culture medium with 2,5 ml of fresh Amniomed Plus.
5. On Day 7 / Day 8 check for progress of growth using an inverted microscope. When the areas of growth show a number of cells in mytosis equal to three to four microscope fields (100X), add 30 µl of Colcemid (10 µg/ml).
6. Incubate at 37°C, 5% CO₂ for 60 - 90 minutes.

English

Chromosome preparation, manual method

1. Suck the medium up completely.
2. Drop 3 ml of hypotonic solution (distilled H₂O, 0,6% sodium citrate, 0,1% KCl) and incubate at room temperature for 10 minutes.
3. Suck the hypotonic solution up and drop 3 ml of Ibraimov solution (distilled H₂O, 5% acetic acid), incubate at room temperature for 5 minutes.
4. Suck the Ibraimov modified solution up and drop 3 ml of fresh fixative mixture (methanol : acetic acid, 4 : 1) at room temperature. The time for this wash has no influence on final results.
5. Repeat 4) twice.
6. Remove the glass slide and dry under constant condition of right temperature and relative humidity into Optichrome (ref.EKAMH950).

Chromosome preparation, automated method

1. Put Amniodish in Autochrome (ref. EKAMH1001).
2. Chose the proper program.
3. At the end of the automated process remove the glass slide and dry under constant condition of temperature and relative humidity into Optichrome.

Chorionic villi sampling (CVS) culture protocol

1. Transfer the specimen into a 60 mm Petri dish containing 5 ml of RPMI 1640.
2. Wash villi with fresh medium to remove blood cells.
3. Using the inverted microscope, carefully dissect any remaining clots or decidual fragments.

In situ culture

1. Transfer the villi in a 15 ml sterile centrifuge tube containing 1 ml of Pronase E (4*10⁶ PU/g) and incubate at room temperature for 4 - 6 minutes.
2. Add 3 - 5 ml of cold Hank's Balanced Salts Solution (kept at +4°C).
3. Centrifuge for 5 minutes at 1500 rpm and discard the supernatant.
4. Add 2 ml of sterile Collagenase type II (1 mg/ml) and incubate at 37°C for 10 minutes.
5. Add 3 - 5 ml of cold HBSS (kept at +4°C).
6. Centrifuge for 5 minutes at 1500 rpm and discard the supernatant.
7. Add 2 - 3 ml of Amniomed Plus and resuspend the pellet.
8. Prepare 6 Amniodish, depending on the number of cells. Add 2 ml of Amniomed Plus medium to each Amniodish.
9. Drop about 500 µl of the cell suspension on each Amniodish.
10. Incubate at 37°C in 5% CO₂.
11. On Day 5 replace the culture medium with 2,5 ml of fresh Amniomed Plus.
12. On Day 6 / Day 7 check for progress of growth using an inverted microscope. When the areas of growth show a number of cells in mytosis equal to three to four microscope fields (100X), add 15 µl of Colcemid (10 µg/ml).
13. Incubate at 37 °C, 5% CO₂ for 4 - 6 hours.

Chromosome preparation, manual or automated method

See protocol chromosome preparation of Amniotic fluid.

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EuroClone S.p.A.