Chromosome *synchro P*

Protocollo di utilizzo
Protocole d'utilisation
Benutzungsprotokoll
Protocol of use
Protocolo de utilización
Protocolo de utilização
Πρωτόκολλο χρήσης

- PER COLTURADI CELLULE
- POUR LA CULTURE DES CELLULES
- FÜR ZELLKULTUREN
- FOR CELL CULTURE
- PARA CULTIVO DE CÉLULAS
- PARA CULTURA DE CÉLULAS
- ΓΙΑ ΚΑΛΛΙΕΡΓΕΙΑ ΚΥΤΤΑΡΩΝ

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Kit for culture and synchronization of peripheral blood cells

Material supplied

EK AMTSY 20 (20 tests)

- Chromosome Synchro P, 20 culture tubes
- 2 vials of synchronization kit - solution A
- 2 vials of synchronization kit - solution B

EK AMTSY 50 (50 tests)

- Chromosome Synchro P, 50 culture tubes
- 4 vials of synchronization kit - solution A
- 4 vials of synchronization kit - solution B

EK AMTBSY 100.2 (nearly 40 tests)

- Chromosome Synchro P, 2 bottles 100 ml
- 4 vials of synchronization kit - solution A
- 4 vials of synchronization kit - solution B

EK AMTBSY 100.5 (nearly 100 tests)

- Chromosome Synchro P, 5 bottles 100 ml
- 8 vials of synchronization kit - solution A
- 8 vials of synchronization kit - solution B

EK AMTBSY 500 (nearly 100 tests)

- Chromosome Synchro P, 1 bottle 500 ml
- 8 vials of synchronization kit - solution A
- 8 vials of synchronization kit - solution B

English

Product stability

Chromosome Synchro P is shipped at 4°C. Once received, the product has to be stored at $-20 \pm 2^\circ\text{C}$. The product shelf life is 12 months. Once thawed, the medium may be stored between 2°C and 8°C for a maximum of 30 days. Repeated freezing / thawing cycles of synchronization kit will not affect performances.

Protocol of use

Culture preparation

1. Thaw Chromosome Synchro P culture tube or bottle at room temperature or at $37 \pm 2^\circ\text{C}$.
2. Aseptically add up to 0.5 ml of whole blood to one test tube or dispense 5 ml of bottled medium in a T-Flask and add the sample.
3. Replace the cap of the culture tube.
4. Mix the culture tube (or the T-Flask) contents and incubate at $37 \pm 2^\circ\text{C}$ placing the tube horizontally (on its flat surface). Addition of CO_2 is not necessary as the medium composition is designed for optimal lymphocyte proliferation in a closed system.

Lymphocyte culture synchronization

1. After 48 - 72 hours of culture, add 0.1 ml of solution A and incubate overnight.
2. After the incubation add 0.1 ml of solution B and incubate for 5 hours (it is not necessary to wash the cells).
3. Add 20 μl of Colcemid (10 $\mu\text{g}/\text{ml}$) and incubate at $37 \pm 2^\circ\text{C}$ for 60 minutes (approx.). This last incubation time may be reduced to 15 minutes if the aim is obtaining a higher number of prometaphases.

Chromosome fixation and slide preparation

1. After incubation, transfer the culture in a tube, if necessary, and centrifuge for 3 - 4 minutes at 2.000 rpm (higher speeds will not damage the cells).
2. Discard the supernatant, taking care not to lose part of the cell pellet (it is recommended to leave approximately 0.5 ml of material in the tube).
3. Vigorously resuspend the pellet.
4. Add 5 ml of hypotonic solution (distilled H_2O , 0,56% KCl) and mix adequately.
5. Incubate for at least 7 minutes at room temperature.
6. Centrifuge at 2.000 rpm for 3 - 4 minutes.
7. Discard the supernatant according to step 2.
8. Vigorously resuspend the pellet.
9. Add 5 ml of Ibraimov solution (distilled H_2O , 5% acetic acid) and mix adequately.
10. Immediately centrifuge at 2.000 rpm for 3 - 4 minutes.
11. Discard the supernatant according to step 2.
12. Vigorously resuspend the pellet.
13. Add 5 ml of fixative solution (methanol : acetic acid, 3 : 1) and mix adequately.
14. Repeat steps 10-13.
15. Immediately centrifuge at 2.000 rpm for 3 - 4 minutes.
16. Carefully discard almost all the supernatant.
17. Add a few drops of freshly prepared fixative, taking into consideration that the number of drops should be appropriate to the cell concentration in the pellet.
18. Drop a small aliquot of the cell suspension onto a clean and wet microscope slide.
19. Evaporate under constant condition of right temperature and relative humidity into Optichrome (Ref. EKAMH950).

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EuroClone S.p.A.
EuroClone S.p.A. has a Quality System certified in compliance