Instructions For Use

REF: LPH029

IGH/MAF Translocation, Dual Fusion

PROFESSIONAL USE ONLY

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
Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

The translocation t(14;16) (q32.3;q22) involving c-MAF and IGH is present in 25% of myeloma cell lines¹ and in 2-6% of primary Multiple Myeloma samples^{2,3} and results in high levels of c-MAF expression. Patients harbouring t(14;16) appear to have a worse outcome.^{4,5} The majority of the breakpoints are dispersed over an approximately 500 kb region centromeric to the c-MAF proto-oncogene at 16q23. However, other breakpoints have been seen telomeric to c-MAF.¹ Translocation with the centromeric breakpoints places c-MAF under control of the strong 3 ϕ IGH enhancer. Recent gene expression profiling of myeloma cell lines revealed that c-MAF transactivated cyclin D2, a promoter of cell cycle progression, thus enhancing proliferation of myeloma cells.⁴ The putative tumour suppressor WWOX gene spans the common chromosomal fragile site 16D (FRA16D) at chromosome 16q23.3-24.1. This region is a frequent target for loss of heterozygosity and chromosomal rearrangements in ovarian, breast, hepatocellular, prostate carcinomas and other neoplasias.⁵ WWOX is contained within the breakpoint region of the c-MAF-IGH translocation.

Probe Specification

The IGH/MAF probe mix consists of a green probes covering the Constant, J, D and Variable segments of the IGH gene, flanking the breakpoint region and red probes flanking the breakpoint region in c-MAF. The probe contig. centromeric to c-MAF covers part of the WWOX gene. In the normal cell, these probes will appear as discreet red and green spots, one for each homologue (resulting in a 2G 2R conformation). In a t(14;16) patient, two (yellow) fusion signals in addition to the red and green signals of the normal chromosomes 16 and 14 respectively (1R 1G 2Y) will appear.

Materials Provided

Probe: 50 μ l per vial (5 tests) or 100 μ l per vial (10 tests)

Amount of IGH probe: 120-150 ng/test

Amount of MAF probe: 40-50 ng/test

The probes are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Counterstain: 150 μ l per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 μ g/ml DAPI (4,6-diamidino-2-phenylindole))

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 μ l - 200 μ l
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x 63 or x 100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three probe fluorophores and DAPI simultaneously.

Sample Preparation

The kits are designed for use on cultured peripheral blood cells and bone marrow fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

Slide Preparation

- Spot cell sample onto cleaned microscope slide.
- Wash slide in 2 x SSC for 2 minutes at room temperature (RT).
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins at RT. Allow to air dry.

Pre-Denaturation

- Remove probe from -20°C freezer and allow to warm to RT.
- Ensure probe solution is uniform by repeated pipette mixing.
- Remove 10 μ l of probe per test, place in a microcentrifuge tube. Return remaining probe to -20°C
- Place probe, sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 10 mins.
- Spot 10 μ l probe mixture onto the cell sample and carefully apply coverslip. Seal with rubber solution glue and allow to dry completely.

Denaturation

- Denature sample and probe simultaneously by heating slide on a hotplate at 75°C (+/- 1°C) for 2 mins.

Hybridisation

- Place slide in humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

- Remove coverslip and all traces of glue carefully.
- Wash slide in 0.4 x SSC (pH 7.0) at 72°C (+/- 1°C) for 2 mins.
- Drain slide and wash in 2 x SSC, 0.05% Tween-20 (pH 7.0) at RT for 30 seconds.
- Drain the slide and apply 10 μ l of DAPI antifade.
- Cover with a coverslip and allow colour to develop in the dark for 10 mins.
- View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail. probes@cytocell.com

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Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
CytoCELL Ltd.