## Instructions For Use

REF: LPH030

## IGH/FGFR3 Translocation, Dual Fusion

PROFESSIONAL USE ONLY

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
Further information available at [www.cytocell.com](http://www.cytocell.com)

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

The t(4;14) (p16.3;q32.3) translocation is found in approximately 10% of myeloma patients and results in the deregulation of at least two genes, MMSET (Multiple Myeloma SET domain containing protein) also known as WHSC1 and FGFR3 (fibroblast growth factor receptor 3) with the formation of a fusion product between MMSET and the immunoglobulin heavy chain (IGH) locus and overexpression of FGFR3. The majority of the breakpoints on chromosome 4 occur within the MMSET locus resulting in a fusion gene where the 5' exons of MMSET are replaced with the VDJ region of the heavy chain locus. FGFR3 lies at least 50 kb from the breakpoints after translocation to chromosome 14 and it has been suggested that FGFR3 expression is switched on by close proximity to the strong enhancer elements in the immunoglobulin heavy chain.<sup>1</sup> Other studies state that the breakpoints on chromosome 4 lie in a 113 kb region between FGFR3 and MMSET exon 5.<sup>2</sup> The breakpoint on chromosome 14 is almost exclusively in or near the switch region of the IGH locus.

### Probe Specification

The IGH/FGFR3 probe mix consists of green probes covering the Constant, J, D and Variable segments of the IGH gene in association with 117 kb red probe covering the region telomeric to FGFR3 and a 125 kb red probe centromeric to MMSET gene on chromosome 4. In the normal cell, these probes will appear as discrete red and green spots, one for each homologue (resulting in a 2G 2R conformation). In a t(4;14) patient, two (yellow) fusion signals in addition to the red and green signals of the normal chromosomes 4 and 14 respectively (1R 1G 2Y) will appear.

### Materials Provided

**Probe:** 50µl per vial (5 tests) or 100 µl per vial (10 tests)

Amount of IGH probe: 120-150 ng/test

Amount of FGFR3 probe: 30-40 ng/test

The probes are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC).

**Counterstain:** 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole))

### Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

### Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

### Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl - 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

### Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x 63 or x 100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three probe fluorophores and DAPI simultaneously.

### Sample Preparation

The kit is designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

### FISH Protocol

#### Slide preparation

1. Spot cell sample onto microscope glass slide. Allow to dry.
2. Wash slide in 2 x SSC for 2 minutes at room temperature (RT).
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins at RT. Allow to dry.

#### Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated gentle, pipette mixing.
6. Remove 10µl of probe per test, place in a microcentrifuge tube. Return remaining probe to -20°C.
7. Place probe, sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 10 mins.
8. Spot 10µl probe mixture onto the cell sample and carefully apply coverslip. Seal with rubber solution glue and allow to dry completely.

#### Denaturation

9. Denature sample and probe simultaneously by heating slide on a hotplate at 75°C (+/- 1°C) for 2 mins.

#### Hybridisation

10. Place slide in humid, lightproof container at 37°C (+/- 1°C) overnight.

#### Post-Hybridisation Washes

11. Remove coverslip and all traces of glue carefully.
12. Wash slide in 0.4 x SSC (pH 7.0) at 72°C (+/- 1°C) for 2 mins.
13. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
14. Drain the slide and apply 10 µl of DAPI antifade.
15. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
16. View with fluorescence microscope.

### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

### Procedural Recommendations

1. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
2. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

### Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail.

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#### Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
**Cytocell Ltd.**