Instructions For Use

REF: LPH034

IGK Breakpart

PROFESSIONAL USE ONLY

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
Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Translocations involving the immunoglobulin loci are recurring events in various subtypes of B cell lymphomas. In addition to translocations involving the IGH locus, variant translocations have been described in 5-10% of B cell neoplasms involving either the immunoglobulin kappa (IGK) light chain locus at 2p11.2 or the lambda light chain (IGL) at 22q11.1. The best known translocations involving IG light chain loci are the variant Burkitt's translocations t(2;8)(p12;q24) and t(8;22)(q24;q11), present in up to 25% of all Burkitt's lymphomas.² Other translocations involve the BCL6 oncogene: the t(2;3)(p12;q27) and t(3;22)(q27;q11) and BCL2 locus: t(2;18)(p12;q21) and t(18;22)(q21;q11).^{3,4} Translocations involving the IG light chain loci usually lead to breakage within the joining region of the respective locus.⁵ IGK has largely duplicated structure with duplicate gene regions being 96-100% identical. The IGK proximal copy is present in a 542 kb contig, with 22 potentially functional variable (IGKV) gene segments and 18 pseudogenes plus five joining (J)-segments and one IGK constant (IGKC) gene segment. The distal copy is a 433 kb contig, with 21 potentially functional IGKV gene segments and 15 pseudogenes. The two contigs. are separated by a DNA sequence of 800 kb without any IGKV segments.

Probe Specification

The IGK product consists of a 344 kb probe set labelled in red spanning the distal IGK variable region and 531 kb green probe spanning part of the proximal IGK variable region as well as the joining segments and a constant segment. In the normal cell, 2 red/green (can appear as yellow) signals are expected, (2Y). A translocation will result in 1 green and 1 red signal from the translocated chromosomes and 1 fused red/green signals (can appear as yellow) from the normal chromosome (1R, 1G, 1Y).

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of IGK red probe: 40-50 ng/test

Amount of IGK green probe: 240-300 ng/test

The probes are provided pre-mixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole))

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 µl - 200 µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0,5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x 63 or x 100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three probe fluorophores and DAPI simultaneously.

Sample Preparation

The kits are designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

- Spot cell sample onto microscope glass slide. Allow to dry.
- Wash slide in 2 x SSC for 2 minutes at room temperature (RT).
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins at RT. Allow to dry.

Pre-Denaturation

- Remove probe from -20°C freezer and allow to warm to RT.
- Ensure probe solution is uniform by repeated gentle, pipette mixing.
- Remove 10µl of probe per test, place in a microcentrifuge tube. Return remaining probe to -20°C.
- Place probe, sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 10 mins.
- Spot 10µl probe mixture onto the cell sample and carefully apply coverslip. Seal with rubber solution glue and allow to dry completely.

Denaturation

- Denature sample and probe simultaneously by heating slide on a hotplate at 75°C (+/- 1°C) for 2 mins.

Hybridisation

- Place slide in humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

- Remove coverslip and all traces of glue carefully.
- Wash slide in 0.4 x SSC (pH 7.0) at 72°C (+/- 1°C) for 2 mins.
- Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
- Drain the slide and apply 10 µl of DAPI antifade.
- Cover with a coverslip and allow colour to develop in the dark for 10 mins.
- View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail. probes@cytocell.com

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Patents and Trademarks

Aquarius and Cytocell are registered trademarks of Cytocell Ltd. Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

Cytocell Ltd.
 4 Technopark
 Newmarket Road
 Cambridge, CB5 8PB, UK.
 T: +44(0)1223 294048
 F: +44(0)1223 294986
 E: probes@cytocell.com
 W: www.cytocell.com