Instructions For Use

REF: LPH036

EV11 Breakpart

PROFESSIONAL USE ONLY

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
Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Chromosome rearrangements involving band 3q26.2 are associated with myeloid malignancies, aberrant expression of the human ecotropic virus integration site-1 (EV11) gene, an unfavourable prognosis and an aggressive clinical course.¹ The EV11 gene encodes a 1051 amino acid finger protein inappropriately expressed in the leukemic cells of 2-5% of AML and MDS patients.² This activation often follows a chromosomal rearrangement involving 3q26 and the two most common translocations: the t(3;3)(q21;q26) and inv(3)(q21q26). EV11 also exists as a longer protein that includes 188 additional amino acids at the N-terminus, named MDS1/EV1. The breakpoints for the translocations and inversions vary considerably. Inv(3) breakpoints are found centromeric to and including the EV11 gene and cover about 600 kb. The breakpoints in 3q26.2 translocations are telomeric to the EV11 gene and cover a region including the telomeric end of the MDS1 gene and the MYNN gene.¹

Probe Specification

Cytocell's product contains a blue probe which covers 562 kb centromeric to the EV11 gene, a green 181 kb probe covering the entire EV11 gene and flanking regions and a red 124 kb probe telomeric to MDS1 and MYNN genes. The three color strategy both shows the presence of a translocation or inversion and allows for each different type of rearrangement to be distinguished. In the normal cell, 2 red/green/blue localised signals are expected. In a breakpoint translocation depending on the breakpoint 1 green/blue signal and 1 separate red from the translocated chromosome and 1 fused red/green/blue signal from the normal chromosome (1R 1B/G and 1R/G/B). Or 1 green/red (yellow) and 1 separate blue signal from the translocated chromosome and 1 fused red/green/blue signal from the normal chromosome (1B 1R/G and 1R/G/B).

Materials Provided

Probe: 50µl per vial (5 tests) or 100 µl per vial (10 tests)
Amount of EV11 red probe: 40-50 ng/test
Amount of EV11 green probe: 120-150 ng/test
Amount of EV11 blue probe: 560-700 ng/test

The probes are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Counterstain:

150 µl per vial (15 tests)
The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole))

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 µl - 200 µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x 63 or x 100. The Triple bandpass filter DAPI/FITC/Texas Red/Aqua is optimal for viewing all three probe fluorophores and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

- Spot cell sample onto microscope glass slide. Allow to dry.
- Wash slide in 2 x SSC for 2 minutes at room temperature (RT).
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins at RT. Allow to dry.

Pre-Denaturation

- Remove probe from -20°C freezer and allow to warm to RT.
- Ensure probe solution is uniform by repeated gentle, pipette mixing.
- Remove 10µl of probe per test, place in a microcentrifuge tube. Return remaining probe to -20°C.

- Place probe, sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 10 mins.

- Spot 10µl probe mixture onto the cell sample and carefully apply coverslip. Seal with rubber solution glue and allow to dry completely.

Denaturation

- Denature sample and probe simultaneously by heating slide on a hotplate at 75°C (+/- 1°C) for 2 mins.

Hybridisation

- Place slide in humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

- Remove coverslip and all traces of glue carefully.
- Wash slide in 0.4 x SSC (pH 7.0) at 72°C (+/- 1°C) for 2 mins.
- Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
- Drain the slide and apply 10 µl of DAPI antifade.
- Cover with a coverslip and allow colour to develop in the dark for 10 mins.
- View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail.
probes@cytocell.com

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.