Instructions For Use

REF: LPH 050

TLX3 Breakpart Probe

PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com


Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

In common with other factors involved in T-ALL, the normal expression of TLX3 (T-Cell Leukaemia Homeobox 3) is disrupted by a chromosomal translocation¹. Unlike TLX1 (HOX11), however, the deregulation of TLX3 is not brought about by close juxtaposition with T-Cell receptor genes, rather, it is brought into contact with another gene which is highly expressed in normal T-cell differentiation, BCL11B (B-Cell Lymphoma 11B, or CTIP2) on chromosome 14q32.1¹. The translocation t(5;14)(q35;q32) is generally cryptic and does not actually disrupt TLX3, rather, it disrupts RANBP17 (Ran-Binding Protein 17) over a 110kb breakpoint region². However, this gene is very close to TLX3 and though RANBP17 expression is not affected by the translocation, TLX3 expression is. Deregulated expression of TLX3 is observed in around 20% of paediatric cases and around 13% of adult T-ALL patients^{3,4}.

Probe Specification

TLX3 region, centromeric 5q35, Red
TLX3 region, telomeric 5q35, Green

The TLX3 probe consists of a red probe (162kb) centromeric to the TLX3 gene and a green probe (173kb) telomeric to the gene. The red probe spans a 162kb proportion of the RANBP17 gene including the D5S653 marker and the green probe spans the NPM1 and FGF18 genes including the marker RH80894. This probe is designed to analyse the occurrence of translocations in this region in T-cell ALL.

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red TLX3 probe: 48-60 ng/test.

Amount of green TLX3 probe: 96-120 ng/test.

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain: 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 24x24 mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights should be limited at all times)

Slide preparation

1. Spot the cell sample onto a glass microscope slide. Allow to dry.
2. Immerse slide in 2xSSC for 2mins at room temperature (RT) without agitation.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2mins at RT.
4. Allow to dry.

Pre-Denaturation

5. Remove the probe from the freezer and allow it to warm to RT.
6. Ensure that the probe solution is uniformly mixed with a pipette.
7. Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to -20°C.
8. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5mins.
9. Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2mins.

Hybridisation

11. Place slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

12. Remove the coverslip and all traces of glue carefully.
13. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2mins without agitation.
14. Drain the slide and immerse in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30secs without agitation.
15. Drain the slide and apply 10µl of DAPI antifade onto each sample.
16. Cover with a coverslip, remove bubbles and allow colour to develop in the dark for 10mins.
17. View with a fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Ltd.

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results

In the normal cell due to the close juxtaposition of the probes expect to see the two homologues as either two fused yellow signals (2Y) or two red/green signals (2RG). In the event of a monoallelic translocation you will see one fused yellow signal (1Y) or a red/green signal (1RG) and one distinct one red (1R) one green (1G) split signal. In the event of a biallelic translocation no fused signals would be present and you will observe two red (2R) and two green (2G) signals.

Customer Support

Please contact the CytoceLL Sales and Marketing Department by telephone or e-mail probes@cytoceLL.com.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.