Instructions For Use

REF: LPE R/G

A range of human alpha and classical satellite probes specific for each chromosome directly labelled in either red or green fluorophore

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Introduction

Satellite probes are specific for human chromosomes; they are highly repeated human DNA sequences for the centromere, pericentromeric or heterochromatic block of each of the 24 chromosomes. The probes will enable the identification and enumeration of human chromosomes in interphase cells or metaphase chromosomes from peripheral blood samples. For detailed probe specifications refer to Table 1.

The probes are produced in a concentrated form to allow mixing, if required, of up to 3 probes in the same hybridisation, from Cytocell's range of concentrated Aquarius satellite probes. A final volume of 10 µl of probe solution is required per hybridisation.

This Aquarius kit contains only one of the probes from the range of directly labelled human alpha and classical satellite probes.

Material Provided

Probe: 15 µl per vial (5 tests)
Amount of satellite probe labelled with green fluorophore: minimum of 20 ng/test
Amount of satellite probe labelled with red fluorophore: minimum of 5 ng/test
The directly labelled probe is provided in hybridisation solution (Formamide; Dextran Sulphate; SSC). It is directly labelled with either a green (FITC spectrum) or a red fluorophore (Texas Red spectrum).

Hybridisation solution (Formamide; Dextran Sulphate; SSC) (150 µl per vial)

Counterstain : (150 µl per vial).

The counterstain is DAPI in antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole))

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formaldehyde which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 µl to 200 µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three probe fluorophores and DAPI simultaneously. For some samples and cells single bandpass filters for FITC, Texas Red or TRITC may be required to optimise signal viewing.

Sample Preparation

The kits are designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures. Baking of slides is not recommended as it may reduce signal fluorescence.

FISH Protocol

Slide Preparation

- Spot cell sample onto cleaned microscope slide.
- Wash slide in 2 x SSC for 2 minutes at room temperature (RT).
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.
- Allow to air dry.

Probe Preparation

- Remove probe and hybridisation solution from 620°C freezer and allow to warm to RT.
- Ensure both solutions are uniform by repeated pipette mixing.
- Using fresh pipette tips remove: final volume of 10 µl of probe solution.
- for a **single probe hybridisation**: 3 µl of probe and 7 µl of hybridisation solution per test
- for a **two probe hybridisation**: 3 µl of each probe and 4 µl of hybridisation solution per test
- for a **three probe hybridisation**: 3 µl of each probe and 1 µl of hybridisation solution per test

And place in a single microcentrifuge tube, gently vortex to mix and pulse spin in a microcentrifuge.

Denaturation

- Prewarm both sample slide and probe at 37°C (+/- 1°C) for 5 minutes.
- Place the 10µl of probe on the spotted area of the slide, carefully apply a 24mm x 24mm coverslip. Seal with rubber solution glue and allow to dry completely.
- Denature at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

- Incubate for **1 hour to overnight** at 37°C (+/- 1°C) in a humid, lightproof container.

Post-Hybridisation Washes

- Remove coverslip and all traces of glue carefully.
- Wash slide in **0.25 x SSC** at 72°C (+/- 1°C) for 2 minutes without agitation*.
- Wash slide in 2 x SSC, 0.05% Tween 20 at RT, pH 7.0 for 30 seconds without agitation.

Counterstain

- Apply 10 µl of the DAPI antifade. Apply a coverslip and allow colour to develop in the dark for 10 minutes. View with a fluorescence microscope.

* If final signal is poor, repeat FISH using 0.4 x SSC post-hybridisation wash.

Stability of Finished Slides

FISHED slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Expected Results for a single probe hybridisation

There will be differences in the relative size of signals observed between chromosomes due to the difference in copy number of repeat sequences between chromosomes. The centromere signals for 4, 16, 19 and 20 are particularly small.

- Using one of the satellite probes for the chromosomes 1-12 and 15-20 (except 1/5/19 probe) a diploid sample should show a fluorescent signal at the centromere of both of the chromosomes of the corresponding chromosome in 70 to 90% of cells analysed.
- Using the probe for chromosomes 1/5/19 a diploid sample should show a fluorescent signal at the centromere of each chromosome for chromosomes 1, 5 and 19 in 70 to 90% of cells analysed.
- Using the probe for chromosomes 13/21 or 14/22 a diploid sample should show a fluorescent signal at the centromere of both of the chromosomes for chromosomes 13 and 21 or 14 and 22 in 70 to 90% of cells analysed.
- The chromosome 1 satellite III probe may show faint cross-hybridisation to the pericentromeric region of chromosome 9. This may be reduced when using a 0.25 x SSC stringent wash, compared to a 0.4 x SSC stringent wash.
- The chromosome 2 α-satellite probe may show faint cross-hybridisation to the centromere of an F group chromosome. This may be reduced when using a 0.25 x SSC stringent wash, compared to a 0.4 x SSC stringent wash.
- The chromosome 4 α-satellite probe may show faint cross-hybridisation to the centromeric region of a C group chromosome. This may be reduced when using a 0.25 x SSC stringent wash, compared to a 0.4 x SSC stringent wash.

Table 1 : Probe Specifications

Chr	Catalogue Number *	Locus	Chromosome Region	DNA Class	Probe brightness **
					Red/ Green
1 _a	LPE 01G/R	D1Z1	1q12	satellite III	3
2 ^b	LPE 02G/R	D2Z2	2p11.1-q11.1	α-satellite	3
3	LPE 03G/R	D3Z1	3p11.1-q11.1	α-satellite	3
4 ^c	LPE 04G/R	D4Z1	4p11-q11	α-satellite	2
1/5/19	LPE 05G/R	D1Z7 DSZ2 D19Z3	1p11.1-q11 5p11.1-q11.1 19p11-q11	α-satellite	3 2 1
6	LPE 06G/R	D6Z1	6p11.1-q11.1	α-satellite	3
7	LPE 07G/R	D7Z1	7p11.1-q11.1	α-satellite	3
8	LPE 08G/R	D8Z2	8p11.1-q11.1	α-satellite	3
9	LPE 09G/R	D9Z3	9q12	satellite III	3
10	LPE 10G/R	D10Z1	10p11.1-q11.1	α-satellite	3
11	LPE 11G/R	D11Z1	11p11.1-q11.1	α-satellite	3
12	LPE 12G/R	D12Z3	12p11.1-q11.1	α-satellite	3
13/21	LPE 13G/R	D13Z1 D21Z1	13p11.1-q11 21p11.1-q11.1	α-satellite	2-3
14/22	LPE 14G/R	D14Z1 D22Z1	14p11.1-q11.1 22p11.1-q11.1	α-satellite	3
15	LPE 15G/R	D15Z4	15p11.1-q11.1	α-satellite	3
16	LPE 16G/R	D16Z2	16p11.1-q11.1	α-satellite	1-2
17	LPE 17G/R	D17Z1	17p11.1-q11.1	α-satellite	3
18	LPE 18G/R	D18Z1	18p11.1-q11.1	α-satellite	3
20	LPE 20G/R	D20Z1	20p11.1-q11.1	α-satellite	1-2
X	LPE 0XG/R	DXZ1	Xp11.1-q11.1	α-satellite	2
Y	LPE 0YcG/R	DYZ3	Yp11.1-q11.1	α-satellite	2-3
Y	LPE 0YqG/R	DYZ1	Yq12	satellite III	2-3

* G specifies a green label, R specifies a red label.

** Probe brightness: 1. Signal visible with triple filter, but may need individual specific filters for accurate analysis. 2. Signal clearly visible with triple filter. 3. Very bright signal with triple filter.

a, b, c: See the "Expected Results" paragraph for further details.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail. probes@cytocell.com

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Patents and Trademarks

Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
 Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

CytoCELL Ltd.
 4 Technopark
 Newmarket Road
 Cambridge, CB5 8PB, UK.
 T: +44(0)1223 294048
 F: +44(0)1223 294986
 E: probes@cytoCELL.com
 W: www.cytoCELL.com