Instructions For Use

REF: LPE NOR

Acro-P-Arm (Acrocentric chromosome p arms)

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Specification

The NOR (Nucleolar Organizer Regions) probe is specific for acrocentric chromosome p arms. It is designed for fluorescence *in situ* hybridisation to metaphase chromosomes from fixed cultured peripheral blood cells.

NORs are a chromosomal region around which the nucleoli develop at the end of mitosis¹. The region contains several tandem copies of ribosomal RNA genes (5.8S, 18S, and 28S) clustered on the short arms of the acrocentric chromosomes 13, 14, 15, 21 and 22^{2,3}.

A NOR can translocate to a terminal region of another chromosome. In rare instances this can be pathogenic, particularly when the translocation has led to a deletion at the tip of recipient chromosome⁴.

The NORs show a high degree of variations in size and staining pattern. If the region is not translated conventional silver staining will not stain the NOR. However, FISH is able to overcome this problem and acrocentric chromosomes can be accurately confirmed⁵.

Materials Provided

Probe: 100 µl per vial (10 tests)

Amount of Acro-P probe labelled with red fluorophore: 20-25 ng/test

The probes are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole))

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl - 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x 63 or x 100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing probe fluorophore and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines.

Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

1. Spot cell sample onto microscope glass slide. Allow to dry.
2. Wash slide in 2 x SSC for 2 minutes at room temperature (RT).

3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins at RT. Allow to dry.
- Pre-Denaturation
4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated gentle, pipette mixing.
6. Remove 10µl of probe per test, place in a microcentrifuge tube. Return remaining probe to -20°C.
7. Place probe, sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 10 mins.
8. Spot 10µl probe mixture onto the cell sample and carefully apply coverslip. Seal with rubber solution glue and allow to dry completely.

Denaturation

9. Denature sample and probe simultaneously by heating slide on a hotplate at 75°C (+/- 1°C) for 2 mins.

Hybridisation

10. Place slide in humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

11. Remove coverslip and all traces of glue carefully.
12. Wash slide in 0.4 x SSC (pH 7.0) at 72°C (+/- 1°C) for 2 mins.
13. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
14. Drain the slide and apply 10 µl of DAPI antifade.
15. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
16. View with fluorescence microscope.

Expected results: up to ten signals which vary in size.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
2. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.