Instructions For Use

REF: LPT R/G

A range of Subtelomere Specific Probes directly labelled in either a red or a green fluorophore

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation, is treated with heat and Formamide to denature the double-stranded DNA, rendering it single-stranded. The target DNA is thus available for annealing to a similarly denatured, single-stranded, fluorescent-labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

The probes are produced in a concentrated form to allow the mixing, if required, of up to 3 probes in the same hybridisation, from Cytocell's range of concentrated Aquarius subtelomere specific probes. A final volume of 10µl of probe solution is required per hybridisation.

This Aquarius kit contains only one of the probes from the range of directly labelled telomere specific probes.

Material Provided

Probe: 15µl per vial (5 tests)
Amount of subtelomere specific probe: minimum of 40 ng/test
The probe is provided in hybridisation solution (Formamide; Dextran Sulphate; SSC). It is directly labelled with either a green (FITC spectrum) or a red fluorophore (Texas Red spectrum).

Hybridisation solution (Formamide; Dextran Sulphate; SSC) (150µl per vial).

Counterstain: (150µl per vial).

The counterstain is DAPI in antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole))

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formaldehyde which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institutions' guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at -20°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1µl to 200µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing both probe fluorophore and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

- Spot cell sample onto cleaned microscope slide
- Wash slide in 2 x SSC for 2 minutes at room temperature (RT).
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins
- Allow to air dry

Probe Preparation

- Remove probe (s) and hybridisation solution from -20°C freezer and allow to warm to RT
- Ensure solutions are uniform by repeated gentle, pipette mixing
- Using fresh pipette tips remove: final volume of 10µl of probe solution
 - for a **single probe hybridisation**: 3µl of probe and 7µl of hybridisation solution per test
 - for a **two probe hybridisation**: 3µl of each probe and 4µl of hybridisation solution per test
 - for a **three probe hybridisation**: 3µl of each probe and 1µl of hybridisation solution per test
- And place in a single microcentrifuge tube, gently vortex to mix and pulse spin in a microcentrifuge

Denaturation

- Prewarm both sample slide and probe to 37°C (+/- 1°C) for 5 minutes
- Place the 10µl of probe on the spotted area of the slide, carefully apply a 24 x 24mm coverslip. Seal with rubber solution glue and allow to dry completely.

- Denature at 75°C (+/- 1°C) for 2 minutes

Hybridisation

- Incubate overnight at 37°C (+/- 1°C) in a humid, lightproof container

Post-Hybridisation Washes:

- Remove coverslip and all traces of glue carefully
- Wash slide in 0.4 x SSC at 72°C (+/- 1°C) pH 7.0 for 2 minutes without agitation
- Wash slide in 2 x SSC, 0.05% tween 20 at RT, pH 7.0 for 30 seconds without agitation

Counterstain

- Apply 10µl of the DAPI antifade. Apply a coverslip, remove bubbles and allow colour to develop in the dark for 10 minutes. View with a fluorescence microscope

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Expected Results for a single probe hybridisation

- A normal sample will show a signal of the telomeric region of each of the relevant chromosome homologues, or 2 signals in interphase cells.
- The chromosome location of the signals is immediately apparent from the cytogenetic characterisation of the chromosomes and thus both balanced and unbalanced abnormalities can be instantly detected and accurately described.
- # The 2q PAC (dj10111017) detects the common polymorphism described by Macina *et al* in 1994. Should a 2q deletion be detected studies have stressed the importance of testing the deletion by FISH using a second probe mapping to the region (Fan *et al* 2001) and to test the parents to determine its inheritance (Knight *et al* 2000). To this end, a second 2q BAC, 172113, D2S447 which is more proximal has been added to the Aquarius Telomere Specific probe range.

Probe Specifications

Probe	Catalogue Number	Clone name	Marker	Accession Numbers	Max. physical distance from telomere (kb)	Cross-hybridisations & polymorphisms
1p	LPT01PR/G	CEB108	CEB108/T7		<300	
1q	LPT01QR/G	160H23	1qtel19	D1S3739	80	
2p	LPT02PR/G	DJ892G20	2ptel27	D2S2983	330	
2q	LPT02QR/G	DJ1011017	2qtel47	D2S2986	240	Polymorphism #⊗
2q NP	LPT02QNP R/G	172113	VII- YRM2112	D2S447	240	Non-polymorphic#
3p	LPT03PR/G	DJ1186B18	3ptel25	D3S4559	450	
3q	LPT03QR/G	196F4	3qtel06	D3S1272	450	
4p	LPT04PR/G	36P21	4ptel04	D4S3360	73	
4q	LPT04QR/G	DJ963k6	4qtel11		275-500	
5p	LPT05PR/G	189N21	5ptel48		Unknown	
5q	LPT05QR/G	240G13	5qtel70	D5S2097	245	
6p	LPT06PR/G	621I1	6ptel48		300	
6q	LPT06QR/G	57H24	6qtel54	D6S2522	280	
7p	LPT07PR/G	109a6		G31341	<255	
7q	LPT07QR/G	2000a5		G31340	<8	
8p	LPT08PR/G	DJ580L5	8ptel91	D8S2333	250	8p with 1p & 3q
8q	LPT08QR/G	489D14	8qtel11	D8S1925	170	
9p	LPT09PR/G	43N6	9ptel30		600	
9q	LPT09QR/G	112N13	9qtel33	D9S2168	65	9q with 10p & 16p, 18p, XqYq
10p	LPT10PR/G	306F7	10ptel45	D10S2488	320	
10q	LPT10QR/G	137E24	10qtel24	D10S2490	270	
11p	LPT11PR/G	DJ908H22	11ptel03	D11S2071	125	11p with 17p
11q	LPT11QR/G	DJ770G7	11qtel38	D11S4974	65	11q with 12q (interstitial)
12p	LPT12PR/G	496A11	12ptel27		Unknown	12p with 6p, 20q
12q	LPT12QR/G	221K18	12qtel87	D12S2343	190	
13q	LPT13QR/G	163C9	13qtel56	D13S1825	170	
14q	LPT14QR/G	DJ820M16	14qtel01	D14S1420	200	
15q	LPT15QR/G	154P1	WI-5214	D15S936	300	
16p	LPT16PR/G	12114	16ptel05	D16S3400	160	
16q	LPT16QR/G	240G10	16qtel48		200	
17p	LPT17PR/G	202L17	17ptel80	D17S2199	60	
17q	LPT17QR/G	362K4	17qtel13	D17S2200	90	17q with 1p, 6q
18p	LPT18PR/G	74G18	VII- YRM2102	D18S552	220	
18q	LPT18QR/G	DJ964M9	18qtel11	D18S1390	290	
19p	LPT19PR/G	DJ546C11	19ptel29		250-500	
19q	LPT19QR/G	F21283	19qtel82		Unknown	
20p	LPT20PR/G	DJ106L1L	20pthy33	D20S502	180	
20q	LPT20QR/G	81F12	20qtel14		50	20q with 6p
21q	LPT21QR/G	63H24	21qtel07	D21S1575	175	
22q	LPT22QR/G	99K24	22qtel31	D22S1726	120	22q with 2q (interstitial)
XpYp **	LPTXYPR/G	839D20		DXYS129	160	
XqYq ***	LPTXYQR/G	225F6		GDB451553/ DXF23751E/ DXYS61	182	
		C8.2/1			88	

See the ⊗Expected Results paragraph for more details.

⊗ See Knight S. and Flint J. (2000) *J Med Genet* 37: 401-409 and Fan Y *et al.* (2001) *Genet in Med* 3 (6) : 416-421 for further information on this polymorphism.

** This probe is specific for the p-arms of both X and Y

*** This probe is specific for q-arms of both X and Y

Customer Support

Please contact the Cytocell Sales and Marketing Department

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Review: Knight and Flint. J Med Genet 2000; 37:401-409

Knight SJL, Lucas S, Regan R, Brenan M, Nicod A, Lawrie NM, Cardy DLN, Nguyen H, Hudson TJ and Flint J (2000). An optimized set of human telomere clones for studying telomere integrity and architecture. Am. J. Hum. Genet. 67 : 320-332

Yao-Shan Fan, Yang Zhang; Marsha Speevak; Sandra Farrell; Jack Jung; and Victoria Siu. (2001) Detection of Submicroscopic aberrations in patients with unexplained mental retardation by fluorescence *in situ* hybridisation using multiple subtelomeric probes. Genetics in Medicine 3 (6) 416-421

Macina, RA; Negorev, DG; Spais, LA; Hu, Xi; Reithman, HC. (1994) Sequence organisation of the human chromosome 2q telomere. Human Molecular Genetics 3 1847-1853

Patents and Trademarks

Aquarius and Cytocell are registered trademarks of Cytocell Ltd.

Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

Cytocell Ltd.
4 Technopark
Newmarket Road
Cambridge, CB5 8PB, UK.
T: +44(0)1223 294048
F: +44(0)1223 294986
E: probes@cytocell.com
W: www.cytocell.com

004/2010-06-28

DS#080/CE