## Instructions For Use

REF: PMP007/ PMP008/ PMP009

### Chromoprobe Multiprobe® - T System

FOR PROFESSIONAL USE ONLY

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Further information available at [www.cytocell.com](http://www.cytocell.com)

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Telomere analysis for chromosome rearrangements is very attractive as the majority of translocations involve chromosome ends and therefore an assay that targets telomeres will identify all of these with high sensitivity regardless of size. With the advent of the first generation of telomere clones (mainly cosmids 35-40 kb) and applying the Multiprobe technology, Cytocell was able to launch Chromoprobe Multiprobe®-T, a FISH based assay which allows every telomere to be investigated using only a single microscope slide per patient. A directly labelled product incorporating a second generation set of telomere PAC/BAC clones and some of the original cosmid clones has now succeeded the first generation product. The strategy involves FISH of the telomere specific probes to metaphase chromosomes. Deletions are detected by the absence of signal from the end of the cognate chromosome and trisomies by the presence of three signals. The chromosome location of the signals is immediately apparent from cytogenetic characterisation of the chromosomes and thus both balanced and unbalanced translocations can be instantly detected and accurately described.

The Multiprobe device is divided into 24 square areas, each containing a telomere specific probe for both the p arm and q arm of a chromosome (except the D and G groups). The probes are reversibly bound to a glass device, a method exclusive to Cytocell, thereby simplifying fluorescence *in situ* hybridisation by removing the need to prepare the probes. Denaturation of probe and target DNA occurs simultaneously under the device once heated. The probes are differentially labelled. The p arm probe is labelled in a fluorophore specific for the FITC spectrum and the q arm probe is labelled in a fluorophore specific for the Texas Red spectrum.

The FISH protocol incorporates co-denaturation and overnight hybridisation, followed by non-formamide rapid washes. Chromoprobe Multiprobe-T provides an effective assay for identifying telomere rearrangements, is highly specific and cost effective.

#### Material Provided

Each kit contains the following components which are sufficient for either 2 (Cat. No. PMP009), 5 (Cat. No. PMP008) or 10 (Cat. No. PMP007) patient samples:

- 2, 5 or 10 Glass Multiprobe devices coated with 42 directly labelled subtelomeric probes (minimum of 20 µg of probe per square).
- 4, 7 or 12 Glass slides printed with a special template
- 500 µl Hybridisation Solution B: Formamide, Dextran Sulphate, SSC
- 500 µl Counterstain Solution : DAPI (ES : 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)), Antifade
- 1 Cytocell Slide Surface Thermometer
- 1 Cytocell Chromoprobe Multiprobe Hybridisation Chamber
- 4, 7 or 12 Large coverslips (24 mm x 60 mm)

#### Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. The hybridisation solution contains formamide which is toxic. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

#### Storage and Handling

The Chromoprobe Multiprobe System kit should be stored at 2-8°C until the expiry date indicated on the kit label. **Do not freeze.** The counterstain vial must be stored in the dark.

#### Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl - 200 µl
- c) Non-stirring water bath with accurate temperature control at 37°C
- d) Water bath with accurate temperature control at 72°C
- e) Microcentrifuge tubes (0.5 ml)
- f) Fluorescence microscope (Please see Optimal microscope and filter set up section)
- g) Plastic or glass coplin jars
- h) Forceps
- i) Fluorescence grade microscope lens immersion oil
- j) Bench top centrifuge

#### Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

#### Sample Preparation

The Chromoprobe Multiprobe System is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative, which should be prepared according to the laboratory or institution guidelines. Prepare blood metaphase spreads on Cytocell Chromoprobe Multiprobe template slides according to Cytocell protocol below. Baking or otherwise ageing of slides is not recommended as it may reduce signal fluorescence.

#### Chromoprobe Multiprobe Protocol

##### 1. Slide Preparation

- i. Clean template slide
- ii. Soak the template slide for 2 minutes in 100% methanol and polish dry with a clean soft tissue. Establish correct cell density  
It is important that the intended sample has a sufficiently high mitotic index to allow detection of chromosome abnormalities. To check the density of the sample, using a micropipette (e.g. a Gilson P10 or P20) pipette 4 µl of the cell suspension onto one of the areas of the spare template slide and allow to air dry. The small volume used means that you usually have to gently touch the slide with the pipette tip to transfer the suspension. Examine by phase contrast microscopy.  
-If the cell density is too high, dilute the suspension with fresh fixative.  
-If the mitotic index is too low, spin down the fixed cell suspension at 160 xg for 10 minutes. Note the volume of supernatant, remove, and re-suspend the cell pellet in a smaller volume of fresh fixative.  
-If cell sample concentration has been altered, spot 4 µl of the concentrated sample onto another square of your test slide and re-examine by phase contrast microscopy.  
**Please Note: 50 µl is the minimum volume required for the Multiprobe device protocol.**
- iii. Quality control of samples  
Samples should be examined for cytoplasm since this will interfere with the final FISH results. One method for reducing cytoplasm is to spot 2 µl of your sample onto the template slide and watch the fixative as it

spreads out. In the normal situation, the fixative will spread to maximum, recede and then evaporate. To clean up any cytoplasm, we have found that effective results are achieved if a fresh drop of 3:1 fixative is allowed to fall onto the spot at the point when the spreading fixative has reached its maximum spread. Allow the second drop of fixative to evaporate and re-examine the spot.

- iv. Spotting of slide  
Pipette 2 µl of cell suspension onto all 24 areas following the method outlined below.

**Please Note:** If samples are of excellent quality, it is important that each spot is allowed to dry before proceeding to the next. If samples exhibit a degree of cytoplasmic contamination, however, then spotting the samples one after the other in the manner described below may help.

The spotting sequence of alternating squares shown below should be adhered to so that the drops of suspension do not interfere with each other as they spread.

1	3	5	7	
10	12	14	16	
17	19	21	X	

Once the first group of drops has air-dried, spot the remaining squares with 2µl drops in the same manner. After the slide has dried, examination of the slide under phase contrast will reveal whether any squares have been missed and they can be respotted.

If the metaphase cells appear overspread, if there are chromosomes randomly scattered, then allow every spot to dry before proceeding to the next.

If upon examination of slide, a square has insufficient cells/metaphases, spotting further drop(s) of suspension can be added to increase the cell density.

##### 2. Preparation of Multiprobe device and template slide

- i. Place the Chromoprobe Multiprobe Hybridisation Chamber in a 37°C water bath and allow to equilibrate to 37°C (+/- 1°C). (This may take up to an hour if the water bath has been switched on from cold).
- ii. Mix the hybridisation solution by repeated pipetting and pre-warm a 30 µl aliquot per Multiprobe device to 37°C (+/- 1°C) (hotplate or incubator). Also pre-warm Multiprobe device to 37°C (hotplate or incubator) by placing device **label side down**.  
**Do not touch the raised boss surfaces of the Multiprobe device.**
- iii. Wash spotted template slide in 2 x SSC for 2 minutes at room temperature (20°C - 25°C).
- iv. Whilst the Multiprobe device is still at 37°C, dehydrate template slide through an ethanol series (2 minutes each at 70%, 85% and absolute ethanol), air dry and place at 37°C (hotplate or incubator) to warm.
- v. Add 1 µl of pre-warmed hybridisation solution (using a P10 micropipette) to each of the raised bosses of the Multiprobe device while it remains at 37°C.

##### 3. Positioning of template slide over Multiprobe

- i. Carefully invert the template slide so that square 1 is located over the top right hand area of the Multiprobe, marked in yellow (Figure 1).
- ii. Make sure that the template slide is carefully aligned over the matching areas on the Multiprobe device. Carefully lower the slide over the Multiprobe so that the drops of hybridisation fluid make contact with the slide. Apply gentle, even pressure to ensure that the hybridisation fluid is spread to the edges of each of the raised areas on the Multiprobe.


Figure 1. Location of telomeric probes on Multiprobe bosses.

- iii. Lift the slide / Multiprobe carefully and invert so that the slide is under the Multiprobe. Make sure the device does not smear across the template slide as this could cause cross-contamination of the probes.
- iv. Place at 37°C (+/- 1°C) (hotplate or incubator) for 10 minutes.

##### 4. Instructions for use of the Cytocell slide surface thermometer

**Please Note:** The temperature of the hotplate should be checked with the Cytocell slide surface thermometer before proceeding to denaturation.

The liquid crystal thermometer must only be used to check the temperature of a hotplate, it must not be used to monitor the hotplate performance over extended time. Place it onto the surface of the hotplate and wait until the different segments stop changing colour. The correct temperature is indicated by a pale green/gold colour. When the segments appear granular and the colours no longer appear uniform and regular, the thermometer should be discarded as it is exhausted. The lifespan of each thermometer should, however, easily be sufficient for a ten Multiprobe kit.

##### 5. Denaturation

A PCR thermal cycler heating block is NOT suitable for use in place of solid bed hotplate for this procedure.

Transfer the slide / Multiprobe device to the 75°C hotplate taking particular care to hold it level (ensure the sample slide is in good contact with the hotplate). Denature on the hotplate at 75°C (+/- 1°C) for **2 minutes**.

##### 6. Hybridisation

Place slide / Multiprobe device in the prewarmed Chromoprobe Multiprobe Hybridisation Chamber, replace the lid and float the chamber in the 37°C (+/- 1°C) water bath (non-stirring) overnight.

**Please note:** Do not seal the lid on the hybridisation chamber.  
Do not place a lid on the water bath.  
Do not hybridise in an incubator.  
Please ensure that the hybridisation chamber is completely dry. (ie. no water or damp tissue inside the chamber).

##### 7. Post hybridisation stringent washes

- i. Preparation of stringent wash solutions  
1. Solution 1: Prepare a coplin / Hellendahl jar containing 0.4 x SSC. Place in a water bath and allow to reach 72°C (+/- 1°C) adjust pH to 7.0.
2. Solution 2 : Prepare a coplin / Hellendahl jar containing 2 x SSC and 0.05% Tween 20. Allow to stand at room temperature (20°C - 25°C).

**Check the temperature and pH of the solutions in coplin jars and adjust if necessary. The pH should be 7.0 when at the correct temperature.**

- ii. Stringent wash steps  
1. Remove the Multiprobe device carefully from the template slide and place the slide in Solution 1 for 2 minutes. (The Multiprobe device cannot be re-used).
2. Place the slide into Solution 2 for 30 seconds.

**Avoid processing more than two Multiprobe slides through the stringency washes at any one time.**

##### 8. Mounting and visualisation of results

- i. Drain off excess solution by placing the slide on its edge on a tissue for around 20 seconds.  
**Do not allow to dry.**
- ii. Add one 20µl drop of DAPI-Antifade solution to both ends of the slide and apply a large (24mm x 60mm) coverslip, as provided, slowly to avoid the formation of air bubbles.
- iii. Blot with a piece of filter paper / tissue.
- iv. Leave for 10 minutes in the dark before viewing by fluorescence microscopy.
- v. Certain types of microscope have slide holders which make it difficult to view the extreme ends of the slide. If this occurs then simply turn the slide through 180° which will help with the viewing of the slide.

The probes used on the Multiprobe device are directly labelled with fluorophores which are light sensitive. Results are improved when the probes are exposed to minimal amounts of light during these procedures; however, it is not necessary to work in the dark.

#### Stability of Finished Slides

FISHED slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

#### Procedural Recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

#### Expected Results

- For all chromosomes except the acrocentric D & G group chromosomes, a normal sample will show 2 green (p-arm) and 2 red (q-arm) signals on the respective ends of the chromosome. For acrocentric chromosomes, a normal sample will show red signals on the q-arms only.
- \* The 2q PAC (d10i1017) detects the common polymorphism described by Macina *et al* in 1994. Should a 2q deletion be detected studies have stressed the importance of testing the deletion by FISH using a second probe mapping to the region (Fan *et al* 2001) and to test the parents to determine its inheritance (Knight *et al* 2000). To this end, the second 2q BAC, 172113, D2S447 which is more proximal and labelled in red, has been incorporated on square Y of the device. It is supplied with a control probe (2pter) labelled in green. It must be stressed, however, that although further studies have been facilitated by the provision of the second clone on square Y of the device, any apparent 2q deletions must be followed-up with testing of other family members to prevent the misinterpretation of a clinically significant distal 2q deletion as a polymorphism or normal variant.

#### Template slide layout with probe cross-hybridisation and polymorphism details.

Template slide square	Probe	Probe label	Cross-hybridisation and polymorphism of probe clones
1	1p	green	
	1q	red	
2	2p	green	
	2q	red	2q polymorphism*⊗
3	3p	green	
	3q	green	
4	4p	green	
	4q	red	
5	5p	green	
	5q	red	
6	6p	green	
	6q	red	
7	7p	green	
	7q	red	
8	8p	green	8p with 1p and 3q
	8q	red	
9	9p	green	
	9q	red	
10	10p	green	9q with 10p and 16p, 18p, XqYq
	10q	red	
11	11p	green	11p with 17p
	11q	red	11q with 12q (interstitial)
12	12p	green	12p with 6p, 20q
	12q	red	
13	13q	red	
14	14q	red	14q with 16p11.2
15	15q	red	
16	16p	green	
	16q	red	
17	17p	green	
	17q	red	17q with 1p, 5q, 6q and 11p
	18p	red	
18	18q	red	
	19p	green	19p occasionally with 20q
19	19q	red	
	20p	green	
20	20q	red	20q with 6p
	21q	red	
21	21q	red	22q with 2q (interstitial)
22	XpYp**	green	
X	XqYq***	red	
	2p	green	
Y	2q NP	red	non-polymorphic*

\* See the "Expected Results" paragraph for more details.

⊗ See Knight S. and Flint J. (2000) *J Med Genet* 37: 401-409 and Fan Y *et al.* (2001) *Genet in Med* 3 (6) : 416- 421 for further information on this polymorphism.

\*\* This probe is specific for the p-arms of both X and Y

\*\*\* This probe is specific for q-arms of both X and Y

#### Customer Support

Please contact the CytoCell Sales and Marketing Department by telephone or e-mail. probes@cytozell.com

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**References/Bibliographie/Bibliographia/Literatur**

1. Knight SJL and Flint J (2000). Perfect endings : a review of subtelomeric probes and their use in clinical diagnosis. *J. Med. Genet.* 37 : 401-409
2. Knight SJL, Lucas S, Regan R, Brennan M, Nicod A, Lawrie NM, Cardy DLN, Nguyen H, Hudson TJ and Flint J (2000). An optimized set of human telomere clones for studying telomere integrity and architecture. *Am. J. Hum. Genet.* 67 : 320-332
3. Knight-Jones E, Knight S, Heussler H, Regan R, Flint J and Martin K (2000). Neurodevelopmental profile of a new dysmorphic syndrome associated with submicroscopic partial deletion of 1p36.3. *Dev Med. & Child Neurology.* 42 : 201-206
4. Yao-Shan Fan; Yang Zhang; Marsha Speevak; Sandra Farrell; Jack Jung; and Victoria Siu. (2001) Detection of Submicroscopic aberrations in patients with unexplained mental retardation by fluorescence *in situ* hybridisation using multiple subtelomeric probes. *Genetics in Medicine* 3 (6) 416-421
5. Macina, RA; Negorev, DG; Spais, LA; Hu, XI; Reithman, HC. (1994) Sequence organisation of the human chromosome 2q telomere. *Human Molecular Genetics* 3 1847-1853

**Patents and Trademarks**

Chromoprobe, Cytocell and Chromoprobe Multiprobe are registered trademarks of Cytocell Ltd. The Chromoprobe principle is covered by international patents WO9314223, EP0623177. The design of the Multiprobe is a registered design, number 2050801 and is also covered by Design Patent No.436,668. Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

**Cytocell Ltd.**  
 4 Technopark  
 Newmarket Road  
 Cambridge, CB5 8PB, UK.  
 T: +44(0)1223 294048  
 F: +44(0)1223 294986  
 E: probes@cytocell.com  
 W: www.cytocell.com

004/2008-03-20

PI016/CE