Instructions For Use

REF: PMP030 /031/032

**Chromoprobe
Multiprobe® - System
ALL v2 Panel**

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.


Chromoprobe Multiprobe ALL Panel

Cytocell's The Chromoprobe Multiprobe ALL Panel has been designed to provide as much information about the patient's genotype as possible using FISH probes applied to a single slide in a single hybridisation experiment. It detects rearrangements that occur primarily in B-cell lineage ALL, though some T-lineage markers have been included. In addition, in common with the CLL panel the panel strategy has been developed to give the maximum amount of information from cells at interphase and in some cases is capable of detecting chromosome rearrangements that are undetectable using standard cytogenetics.

FISH Tests on the Chromoprobe Multiprobe ALL Panel


MYC

For the MYC fission product, the probes are spaced widely apart (717 kb) so as to be capable of picking up each of the three translocations. The proximal probe (177 kb) is 332 kb 3' of the gene and is labelled in red, whilst the distal probe (188 kb) is 380 kb 5' of the gene and is labelled in green.


P16

The P16 probe, labelled in red, covers a 101 kb region of 9p21, extending from 59 kb 3' of P16 to the 5' end of P15. The probe mix also contains a control probe for chromosome 9 (D9Z3, the heterochromatic block at 9q12), labelled in green.


E2A

The E2A probe mix consists of a 3' probe, labelled in green, covering 164 kb 3' of the gene spanning the markers D19S883 and D2S2123 whilst the 5' region probe, labelled in red, (191 kb) is 157 kb from the 5' end of the gene.



MLL

The MLL probe mix has a green 5' probe and a red 3' probe. The 5' region probe (96 kb) spans the CDG3 gene and 16 kb of the UBE4A gene whilst the 3' region probe (87 kb) spans a region 3' of exon 5.



TEL (ETV6)/AML1

The ETV6 probe mix contains a probe 5' of ETV6 covering a region between the marker D12S845 and the 5' end of intron 2, measuring 182 kb and a second probe covering a region 3' of the gene extending 170 kb from the marker D12S1898. Both are labelled in red. For AML1, the probes cover 151 kb 3' of RUNX1, including CLIC6 and a second probe extending from intron 3 of RUNX1 to 51 Kb beyond the 5' end of the gene. In the normal cell, these probes will appear as discrete red and green spots, one for each homologue (resulting in a 2G 2R conformation). This probe set will detect the t(12;21) fusion and may also detect larger deletions of the TEL1 (ETV6) gene.


BCR/ABL

The BCR probe mix contains a probe 3' of BCR covering a region extending 171 kb 3', covering the genes GNAZ and RAB36 and a second probe covering a region 262 kb 5' of the gene extending 148 kb. Both are labelled in red and are orientated such that breakpoints in either mBCR or MBCR will result in a fusion signal. For ABL, the probe contig covers 349 kb from the middle of the FUB3 gene to a point 64 kb from the 5' end of ABL.

CHIC2


The CHIC2 probe has been provided on the panel as an enumerator for chromosome 4 to assess the levels of hyperdiploidy in ALL patients.

IGH

With all these rearrangements having breakpoints within the IGH gene, we have designed a split probe set for IGH, which can detect any rearrangement, which involves the splitting of the IGH gene in the region between the Constant and Variable segments, thereby identifying the IGH translocation partner chromosome in the less common rearrangements of this gene.


The Constant region of the gene is labelled in red whilst the Variable segment is labelled in green. The normal situation is represented by fusion or close-juxtaposition of the red and green signals (2Y) whilst a rearrangement of the gene is detectable by separate green and red signals (1Y, 1G, 1R). This probe set will detect such rearrangements in interphase cells as well as in dividing cells.

Hyperdiploidy

The ALL panel has been designed so that hyperdiploidy can be determined for chromosomes 4, 10 and 17 using CHIC2 (see above) for chromosome 4 along with centromeric probes for chromosomes 10 and 17. In addition, probes for IGH and AML1 will enumerate chromosomes 14 and 21 (squares 8 and 5 respectively), whilst MYC on square 1 will enumerate chromosome 8 which is associated with a poor outcome in adult ALL. Enumeration of chromosomes 4, and 10 using this panel will identify those patients with a relatively better outcome, as these have been associated with the subgroup of patients exhibiting between 55 and 67 chromosomes.

Please note Chromosome 10 is labelled in Green and Chromosome 17 is labelled in blue.

Material Provided

Each kit contains the following reagents, which are sufficient for either 2 (Cat. No. PMP030), 5 (Cat. No. PMP031) or 10 (Cat. No. PMP032) patient samples:

- 2, 5 or 10 Chromoprobe Multiprobe – Leukaemia Panel devices coated with directly labelled single copy probes
- 4, 7 or 12 Glass slides printed with a special template
- 500 µl Hybridisation Solution B: Formamide, Dextran Sulphate, SSC
- 500 µl Counterstain Solution: DAPI (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)), Antifade
- 1 Cytocell Slide Surface Thermometer
- 1 Cytocell Chromoprobe Multiprobe Hybridisation Chamber

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. The hybridisation solution contains Formamide, which is toxic. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Chromoprobe Multiprobe System kit should be stored at 2-8°C until the expiry date indicated on the kit label. Do not freeze. The counterstain vial must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate with accurate temperature control up to 80°C
- b) 37°C incubator
- c) Variable volume micropipettes range 1 µl -200 µl
- d) 37°C water bath (without stirrer)
- e) Microcentrifuge tubes (0.5 ml)
- f) Water bath with accurate temperature control at 72°C
- g) Fluorescence microscope (Please see Optimal microscope and filter set up)
- h) Plastic or glass coplin jars
- i) Centrifuge
- j) Forceps
- k) Fluorescence grade microscope lens immersion oil
- l) Fluorescence grade glass coverslips (24 x 50 mm)

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FTTC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The Chromoprobe Multiprobe System is designed for use on cultured peripheral blood or bone marrow cells fixed in Carnoy's fixative, which should be prepared according to the laboratory or institution guidelines.


Prepare samples on Cytocell Chromoprobe Multiprobe template slides according to Cytocell protocol below.

Baking or otherwise ageing slides is not recommended as it may reduce signal fluorescence.

Chromoprobe Multiprobe Protocol

I. Slide preparation

- i. Clean template slide
- ii. Soak the template slide for 2 minutes in 100% methanol and polish dry with a clean soft tissue.
- iii. Establish correct mitotic index
It is important that the intended sample has a sufficiently high mitotic index to allow detection of chromosome abnormalities. To check the density of the sample, using a micropipette (e.g. a Gilson P10 or P20) pipette 4 µl of the cell suspension onto one of the areas of the spare template slide and allow to air dry. The small volume used means that you usually have to gently touch the slide with the pipette tip to transfer the suspension. Examine by phase contrast microscopy. If the cell density is too high, dilute the suspension with fresh fixative. If the mitotic index is too low, spin down the fixed cell suspension at 160 xg for 10 minutes. Note the volume of supernatant, remove, and re-suspend the cell pellet in a smaller volume of fresh fixative. If cell sample concentration has been altered, spot 4 µl of the concentrated sample onto another square of your test slide and re-examine by phase contrast microscopy.
Please Note: 50 µl is the minimum volume required for the protocol.
- iii. Quality control of samples
Samples should be examined for cytoplasm since this will interfere with the *in situ* protocol. If the chromosomes appear to be enclosed by a granular material when examined under phase contrast microscopy, then this will compromise results. One method for reducing cytoplasm is to spot 4 µl of your sample onto the template slide and watch the fixative as it spreads out. In the normal situation, the fixative will spread to maximum, recede and then evaporate. To clean up any cytoplasm we have found that effective results are achieved if a fresh drop of 3:1 fixative is allowed to fall onto the spot at the point when the spreading fixative has reached its maximum. Allow the drop of fixative to evaporate and re-examine the spot.
- iv. Spotting of slide
Pipette 4 µl of cell suspension onto all 8 areas of the template slide in a sequence of alternating squares as shown below. This will prevent the cell spreads from interfering with each other.

Once the first group of drops has air-dried, spot the remaining squares with 4 µl drops in the same manner. After the slide has dried, examination of the slide under phase contrast will reveal whether any squares have been missed. If spots have been missed, or squares have too few cells, simply spot those squares again: it is not necessary to re-spot a new slide.

If upon examination of slide, a square has insufficient cells/metaphases, further drop(s) of suspension can be added to increase the cell density.

Please note: If the metaphase cells appear overspread then clean the template slide thoroughly in methanol and re-spot allowing every spot to dry before proceeding to the next.

2. Preparation of ALL Panel device and template slide

- i. Ensure that the Chromoprobe Multiprobe Hybridisation Chamber is in the 37°C water bath and allow equilibrating to 37°C (+/- 1°C), this may take up to an hour if the water bath has been switched on from cold
- ii. Mix the hybridisation solution by repeated pipetting and pre-warm a 25 µl aliquot per device to 37°C. Also pre-warm each device to 37°C by placing the device label side down.
Do not touch the raised boss surfaces of the device.
- iii. Wash template slides containing fixed samples in 2 x SSC for 2 minutes at room temperature (20 - 25°C).
- iv. Whilst the device is still at 37°C, dehydrate template slides containing fixed samples through an ethanol series (2 minutes each in 70%, 85% and absolute ethanol), dry and place at 37°C to warm.
- v. Add 2 µl of pre-warmed hybridisation solution to each of the eight areas on the pre-warmed device using a P10 micropipette while it remains at 37°C.

3. Positioning of template slide over the device

- i. Carefully invert the template slide over the device such that the number 1, which is now upside down, is located over the top right hand area of the device (Figure 1).

To help locate square 1, its position on the device has been marked with a coloured label.


Figure 1. Location of probes on the ALL Panel

- ii. Make sure that the template slide is carefully aligned with the matching areas on the device. Carefully lower the slide over the device so that the drops of hybridisation solution make contact with the slide. Apply gentle, even pressure to ensure that the hybridisation solution is spread to the edges of each of the raised areas on the device.
- iii. Lift the slide carefully holding the frosted end of the glass slide and invert so that the slide is underneath the device. Make sure the device does not smear across the template slide as this could cause cross-contamination of the probes.
- iv. Place at 37°C (+/- 1°C) (hotplate or incubator) for 10 minutes.

4. Instructions for use of the Cytocell Slide Surface Thermometer

The temperature of the hotplate should be checked with the Cytocell Slide Surface Thermometer before proceeding to denaturation.

This thermometer is a liquid crystal device and although reversible, it must be treated with care to ensure a reasonable life span. The thermometer must only be used to check the temperature of a hotplate; it must not be used to monitor the hotplate performance over time. To use the thermometer properly, place it onto the surface of the hotplate and wait until the different segments stop changing colour. The correct temperature is indicated by a pale green / gold colour. When the segments appear granular and the colours no longer appear uniform and regular, the thermometer should be discarded as it is exhausted. The life span of each thermometer should, however, easily be sufficient for a ten-device kit.

5. Denaturation

A PCR thermal cycler-heating block is NOT suitable for use in place of solid bed hotplate for this procedure. Transfer the slide/ device to the hotplate taking particular care to hold it level. Ensure the sample slide is in good contact with the hotplate. Denature on the hotplate at 75°C (+/- 1°C) for 2 minutes.

6. Hybridisation

Place slide/ device in the Chromoprobe Multiprobe Hybridisation Chamber supplied, replace the lid and float the chamber in the 37°C (+/- 1°C) water bath (non - stirring) overnight.

Please note: Do not seal the lid on the hybridisation chamber.

- Do not place a lid on the water bath.
- Do not hybridise in an incubator.
- Please ensure that the hybridisation chamber is completely dry. (ie. no water or damp tissue inside the chamber).

The humidity inside the chamber is vital for optimal hybridisation. The correct levels will be achieved following those steps.

7. Post-hybridisation stringent washes

- i. Preparation of stringent wash solutions
 1. Solution 1: Prepare a Coplin/Hellendahl jar containing 0.4 x SSC. Place in a water bath and allow to reach 72°C (+/- 1°C) adjust pH to 7.0.
 2. Solution 2: Prepare a Coplin/Hellendahl jar containing 2 x SSC and 0.05% Tween 20. Allow to stand at room temperature.

Check the temperature and pH of the solutions in Coplin/Hellendahl jars and adjust if necessary. The pH should be 7.0 when at the correct temperature.
- ii. Stringent wash steps
 1. Remove the device carefully from the slide and place the slide in Solution 1 for 2 minutes. The device cannot be re-used.
 2. Place the slide into Solution 2 for 30 seconds.

Avoid processing more than two slides through the stringency washes at any one time.

8. Mounting and visualisation of results

- i. Counterstaining with DAPI
 1. Apply 20 µl of DAPI to each end of the slide and apply a coverslip (24 x 50 mm)
 2. Blot the slide with filter paper or tissue.
 3. Leave in the dark for 10 minutes before viewing by fluorescence microscopy.
- ii. Certain types of microscope have slide holders, which make it difficult to view the extreme ends of the slide. If this occurs then simply turn the slide through 180°, which will help with the viewing of the slide.

The probes used on the Multiprobe device are directly labelled with fluorophores, which are light sensitive. Results are improved when the probes are exposed to minimal amounts of light during these procedures; however, it is not necessary to work in the dark.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking or otherwise ageing of slides is not recommended as it may reduce signal fluorescence
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, water baths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Expected Results

MYC The normal situation is represented by fusion or close-juxtaposition of the red and green signals (2Y) whilst a rearrangement of the gene is detectable by separate green and red signals (1Y, 1G, 1R). This probe set will detect rearrangements of MYC in interphase cells as well as in dividing cells.

P16 In the normal cell, there should be two red and two green signals (2R, 2G) whilst a deleted cell has either one red signal and two green controls (1R, 2G) if the deletion is hemizygous or no red signals and two green (0R, 2G) if the deletion is homozygous.

E2A. E2A is represented by a break apart (fission) probe. In a normal cell, therefore, there will be two fusion signals (yellow). If E2A is rearranged by translocation, one yellow probe will be split into one red and one green representing the rearranged 19 and the other partner chromosome. The signal pattern will be 1Y, 1R, and 1G.

TEL/AML1. In a t(12;21) patient, there should be two (yellow) fusion signals in addition to the red and green signals of the normal chromosomes 12 and 21 respectively (1R 1G 2Y). Where there are large deletions of ETV6, the signal pattern will be 1R, 2G in the absence of the t(12;21) or 1Y, 1G in the presence of the t(12;21) rearrangement.

BCR/ABL. In the normal cell, these probes will appear as discrete red and green spots, one for each homologous (resulting in a 2G 2R conformation). In a t(9;22) patient, there should be two (yellow) fusion signals in addition to the red and green signals of the normal chromosomes 22 and 9 respectively (1R 1G 2Y).

IGH. IGH is represented by a break apart (fission) probe. In a normal cell, therefore, there will be two fusion signals (yellow). If IGH is rearranged by translocation, one yellow probe will be split into one red and one green representing the rearranged 14 and the other partner chromosome. The signal pattern will be 1Y, 1R, and 1G.

MLL. MLL is also represented by a break apart (fission) probe. In a normal cell, therefore, there will be two fusion signals (yellow). If MLL is rearranged by translocation, one yellow probe will be split into one red and one green representing the rearranged 11 and the other partner chromosome. The signal pattern will be 1Y, 1R, and 1G.

Hybridized ploidy. Signals are visible on a triple filter but for optimal results an Aqua single filter should be used.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail. probes@cytocell.com

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Patents and Trademarks

Chromoprobe, CytoCELL and Chromoprobe Multiprobe are registered trademarks of CytoCELL Ltd.

The Chromoprobe principle is covered by international patents WO9314223, EP0623177. The design of the Multiprobe is a registered design and is also covered by a Design Patent No. 420,745.

Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

CytoCELL Ltd.
4 Technopark
Newmarket Road
Cambridge, CB5 8PB, UK.
T: +44(0)1223 294048
F: +44(0)1223 294986
E: probes@cytoCELL.com
W: www.cytoCELL.com

002/2008-03-20

PI029/CE