Instructions For Use

REF: PMP025/026/027

**Chromoprobe
Multiprobe®- System
AML/MDS Panel**

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
Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

FISH Tests on the Chromoprobe Multiprobe AML/MDS Panel


del(5q)

The EGR1 probe, labelled in red, covers a 167 kb region of 5q31.1, extending from 91.9 kb 3' of EGR1 to 68.4 kb 5' of the gene. The probe mix also contains a control probe for chromosome 5 at 5p15.31 flanking the marker DSS630. In the normal cell, there should be two red and two green signals (2R, 2G). Deleted cells have either one red signal and two green controls (1R, 2G) if the deletion is hemizygous or no red signals and two green (0R, 2G) if the deletion is homozygous.


PML/RAR α

The PML probe mix contains a probe 107 kb proximal to PML, measuring 152 kb, and a second probe 169 kb distal to PML, measuring 107 kb. Both are labelled in red. For RAR α , the probes cover a 194 kb region proximal to RAR α , including CAB3 and a second probe including 13Kb of the distal end of RAR α and extending 153 kb beyond the end of the gene. Both are labelled in green. In the normal cell, these probes will appear as discrete red and green spots, one for each homologue (resulting in a 2G 2R conformation). This probe set will detect the t(15;17) fusion. In a t(15;17) patient, there should be two (yellow) fusion signals in addition to the green and red signals of the normal chromosomes 15 and 17 respectively (1R 1G 2Y).



Tp53

The TP53 probe is 83 kb, labelled in red, and covers a region including 15 kb of the 5' end of TP53 extending 67 kb proximal to the gene, to just beyond the marker D175655. The probe mix also contains a control probe for the 17 centromere (D17Z1). In the normal cell, there should be two red and two green signals (2R, 2G) whilst a deleted cell has either one red signal and two green controls (1R, 2G) if the deletion is hemizygous, or no red signals and two green (0R, 2G) if the deletion is homozygous.

AML1/ETO

The AML1 probe mix contains a probe 151 kb 3' of RUNX1, including CLIC6 and a second probe extending from intron 3 of RUNX1 to 51 kb beyond the 5' end of the gene. Both are labelled in red. For ETO, the probes cover a region 151 kb 5' of ETO, including the marker D8S1950 and a second probe extending from intron 7 of ETO to 127 kb beyond the 3' end of the gene. In the normal cell, these probes will appear as discrete red and green spots, one for each homologue (resulting in a 2G 2R conformation). This probe set will detect the t(8;21) fusion. In a t(8;21) patient, there should be two (yellow) fusion signals in addition to the green and red signals of the normal chromosomes 8 and 21 respectively (1R 1G 2Y).


TRISOMY 8

The AML/MDS panel has been designed so that trisomy 8 can be determined using ETO (see above) on square 4.

MLL


The MLL probe mix also has a green 5' probe and a red 3' probe. The 5' region probe spans the CD3G gene and 16 kb of the UBE4A gene whilst the 3' region probe spans a region 5' of exon 5. The normal situation is represented by fusion or close-juxtaposition of the red and green signals (2Y) whilst a rearrangement of the gene is detectable by separate green and red signals (1Y, 1G, 1R). This probe set will detect rearrangements of MLL in interphase cells as well as in dividing cells.

-7/del(7q)


The proximal 7q probe, labelled in red, covers a 368 kb region of 7q22.1, including 235kb of the 5' end of RELN to beyond the marker D7S658. The distal 7q probe, labelled in green, covers a 203 kb region of 7q31 encompassing TES. In the normal cell, there should be two red and two green signals (2R, 2G). Deleted cells may show one of the following signal patterns:-

- 1 - one red signal and two green (1R, 2G) if the deletion encompasses the proximal CDR only and is hemizygous or
- 2 - no red signals and two green (0R, 2G) if the deletion encompasses the proximal CDR only and the deletion is homozygous or
- 3 - one green signal and two red (1G, 2R) if the deletion encompasses the distal CDR only and is hemizygous or
- 4 - no green signals and two red (0G, 2R) if the deletion encompasses the distal CDR only and the deletion is homozygous or
- 5 - one red signal and one green (1R, 1G) for monosomy 7 or hemizygous deletion of both CDRs on 7q.

CBFB/MYH11


The CBFB probe, labelled in red, covers a 617 kb region of 16q22, extending from 298 kb 5' of CBFB to 246 kb beyond the 3' terminus of the gene. For MYH11, the probe covers 610 kb region of 16p13. In the normal cell, these probes will appear as discrete red and green spots, one for each homologue (resulting in a 2G 2R conformation). This probe set will detect the inv(16) fusion and the variants. In a inv(16) patient, there should be two (yellow) fusion signals in addition to the green and red signals of the normal chromosome 16 (1R 1G 2Y).

Del(20q) (MYBL2)

The proximal 20q probe, labelled in red, covers a 301 kb region of 20q12, extending from 12.4kb inside PTPRT. For the MYBL2 gene in band 20q13.12 the probe contig covers 321kb from 60kb 5' of MYBL2 to 211kb 3' of the gene. In the normal cell, there should be two red and two green signals (2R, 2G). Deleted cells may show one of the following signal patterns:-

- 1 - one red signal and two green (1R, 2G) if the deletion is interstitial encompassing band q12 but not q13.12 and is hemizygous or
- 2 - no red signals and two green (0R, 2G) if the deletion is interstitial encompassing band q12 but not q13.12 and is homozygous or
- 3 - one green signal and two red (1G, 2R) if the deletion encompasses band q13.12 but not q12 and is hemizygous or
- 4 - no green signals and two red (0G, 2R) if the deletion encompasses band q13.12 but not q12 and is homozygous or
- 5 - one red signal and one green (1R, 1G) for monosomy 20 or hemizygous deletion of both bands on 20q.

Material Provided

Each kit contains the following reagents, which are sufficient for either 2 (Cat. No. PMP025), 5 (Cat. No. PMP026) or 10 (Cat. No. PMP027) patient samples:

- 2, 5 or 10 Chromoprobe Multiprobe – Leukaemia Panel devices coated with directly labelled single copy probes.
- 4, 7 or 12 Glass slides printed with a special template
- 500 µl Hybridisation Solution B: Formamide, Dextran Sulphate, SSC
- 500 µl Counterstain Solution: DAPI (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)), Antifade
- 1 Cytocell Slide Surface Thermometer
- 1 Cytocell Chromoprobe Multiprobe Hybridisation Chamber

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. The hybridisation solution contains Formamide, which is toxic. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Chromoprobe Multiprobe System kit should be stored at 2-8°C until the expiry date indicated on the kit label. Do not freeze. The counterstain vial must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate with accurate temperature control up to 80°C
- b) 37°C incubator
- c) Variable volume micropipettes range 1µl - 200µl
- d) 37°C water bath (without stirrer)
- e) Microcentrifuge tubes (0.5 ml)
- f) Water bath with accurate temperature control at 72°C
- g) Fluorescence microscope
- h) Plastic or glass coplin jars
- i) Centrifuge
- j) Forceps
- k) Fluorescence grade microscope lens immersion oil
- l) Fluorescence grade glass coverslips (24 x 50 mm)

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The Chromoprobe Multiprobe System is designed for use on cultured peripheral blood or bone marrow cells fixed in Carnoy's fixative, which should be prepared according to the laboratory or institution guidelines. Prepare samples on Cytocell Chromoprobe Multiprobe template slides according to Cytocell protocol below. Baking or otherwise ageing slides is not recommended as it may reduce signal fluorescence.

Chromoprobe Multiprobe Protocol


1. Slide preparation

- i. Clean template slide
Soak the template slide for 2 minutes in 100% methanol and polish dry with a clean soft tissue.
- ii. Establish correct mitotic index
It is important that the intended sample has a sufficiently high mitotic index to allow detection of chromosome abnormalities. To check the density of the sample, using a micropipette (e.g. a Gilson P10 or P20) pipette 4 µl of the cell suspension onto one of the areas of the spare template slide and allow to air dry. The small volume used means that you usually have to gently touch the slide with the pipette tip to transfer the suspension. Examine by phase contrast microscopy.
If the cell density is too high, dilute the suspension with fresh fixative.
If the mitotic index is too low, spin down the fixed cell suspension at 160 xg for 10 minutes. Note the volume of supernatant, remove, and re-suspend the cell pellet in a smaller volume of fresh fixative.
If cell sample concentration has been altered, spot 4µl of the concentrated sample onto another square of your test slide and re-examine by phase contrast microscopy.
Please Note: 50 µl is the minimum volume required for the protocol.
- iii. Quality control of samples

Samples should be examined for cytoplasm since this will interfere with the *in situ* protocol. If the chromosomes appear to be enclosed by a granular material when examined under phase contrast microscopy, then this will compromise results. One method for reducing cytoplasm is to spot 4 µl of your sample onto the template slide and watch the fixative as it spreads out. In the normal situation, the fixative will spread to maximum, recede and then evaporate. To clean up any cytoplasm we have found that effective results are achieved if a fresh drop of 3:1 fixative is allowed to fall onto the spot at the point when the spreading fixative has reached its maximum. Allow the drop of fixative to evaporate and re-examine the spot.

iv. Spotting of slides

Pipette 4 µl of cell suspension onto all 8 areas of the template slide in a sequence of alternating squares as shown below. This will prevent the cell spreads from interfering with each other.

Once the first group of drops has air-dried, spot the remaining squares with 4µl drops in the same manner. After the slide has dried, examination of the slide under phase contrast will reveal whether any squares have been missed.

If spots have been missed, or squares have too few cells, simply spot those squares again: it is not necessary to re-spot a new slide.

If upon examination of slide, a square has insufficient cells/metaphases, further drop(s) of suspension can be added to increase the cell density.

Please note: If the metaphase cells appear overspread then clean the template slide thoroughly in methanol and re-spot allowing every spot to dry before proceeding to the next.

2. Preparation of AML/MDS Panel device and template slide

- i. Ensure that the Chromoprobe Multiprobe Hybridisation Chamber is in the 37°C water bath and allow equilibrating to 37°C (+/- 1°C). (This may take up to an hour if the water bath has been switched on from cold).
- ii. Mix the hybridisation solution by repeated pipetting and pre-warm a 25 µl aliquot per device to 37°C. Also pre-warm each device to 37°C by placing the device label side down on a 37°C hotplate.
Do not touch the raised boss surfaces of the device.
- iii. Wash template slides containing fixed samples in 2 x SSC for 2 minutes at room temperature (20 - 25°C).
- iv. Whilst the device is still at 37°C, dehydrated template slides containing fixed samples through an ethanol series (2 minutes each in 70%, 85% and absolute ethanol), dry and place at 37°C to warm.
- v. Add 2 µl of pre-warmed hybridisation solution to each of the eight areas on the pre-warmed device using a P10 micropipette while it remains at 37°C.

3. Positioning of template slide over the device

- i. Carefully invert the template slide over the device such that the number 1, which is now upside down, is located over the top right hand area of the device (Figure 1).

To help locate square 1, its position on the device has been marked with a green coloured label.


Figure 1. Location of probes on the AML/MDS Panel

- ii. Make sure that the template slide is carefully aligned with the matching areas on the device. Carefully lower the slide over the device so that the drops of hybridisation solution make contact with the slide. Apply gentle, even pressure to ensure that the hybridisation solution is spread to the edges of each of the raised areas on the device.
- iii. Lift the slide carefully holding the frosted end of the glass slide and invert so that the slide is underneath the device. Make sure the device does not smear across the template slide as this could cause cross-contamination of the probes.
- iv. Place at 37°C (+/- 1°C) (hotplate or incubator) for 10 minutes.

4. Instructions for use of the Cytocell Slide Surface Thermometer

The temperature of the 75°C hotplate should be checked with the Cytocell Slide Surface Thermometer before proceeding to denaturation.

This thermometer is a liquid crystal device and although reversible, it must be treated with care to ensure a reasonable life span. The thermometer must only be used to check the temperature of a hotplate; it must not be used to monitor the hotplate performance over time.

To use the thermometer properly, place it onto the surface of the hotplate and wait until the different segments stop changing colour. The correct temperature is indicated by a pale green / gold colour. When the segments appear granular and the colours no longer appear uniform and regular, the thermometer should be discarded as it is exhausted. The life span of each thermometer should, however, easily be sufficient for a ten-device kit.

5. Denaturation

A PCR thermal cycler-heating block is NOT suitable for use in place of solid bed hotplate for this procedure. Transfer the slide/device to the hotplate taking particular care to hold it level. (Ensure the sample slide is in good contact with the hotplate). Denature on the hotplate at 75°C (+/- 1°C) for 2 minutes.

6. Hybridisation

Place slide/device in the Chromoprobe Multiprobe Hybridisation Chamber supplied, replace the lid and float the chamber in the 37°C (+/- 1°C) water bath (non - stirring) overnight.

Please note: Do not seal the lid on the hybridisation chamber.

Do not place a lid on the water bath.

Do not hybridise in an incubator.

Please ensure that the hybridisation chamber is completely dry. (i.e. no water or damp tissue inside the chamber).

The humidity inside the chamber is vital for optimal hybridisation. The correct levels will be achieved following those steps.

7. Post-hybridisation stringent washes

- i. Preparation of stringent wash solutions
 1. Solution 1: Prepare a Coplin/Hellendahl jar containing 0.4 x SSC. Place in a water bath and allow to reach 72°C (+/- 1°C) adjust pH to 7.0.
 2. Solution 2: Prepare a Coplin/Hellendahl jar containing 2 x SSC and 0.05% Tween 20. Allow to stand at room temperature.

Check the temperature and pH of the solutions in Coplin/Hellendahl jars and adjust if necessary. The pH should be 7.0 when at the correct temperature.

ii. Stringent wash steps

1. Remove the device carefully from the slide and place the slide in Solution 1 for 2 minutes. (The device cannot be re-used)
2. Place the slide into Solution 2 for 30 seconds. Drain off the excess liquid but do not allow to dry.
Avoid processing more than two slides through the stringency washes at any one time.

8. Mounting and visualisation of results.

- i. Counterstaining with DAPI
 1. Apply 20 µl of DAPI to each end of the slide and apply a coverslip (24 x 50 mm)
 2. Blot the slide with filter paper or tissue.
 3. Leave in the dark for 10 minutes before viewing by fluorescence microscopy.
- ii. Certain types of microscope have slide holders, which make it difficult to view the extreme ends of the slide. If this occurs then simply turn the slide through 180°, which will help with the viewing of the slide.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, water baths, and incubators, as these temperatures are critical for optimum product performance.
2. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department by telephone or e-mail. probes@cytocell.com

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Patents and Trademarks

Chromoprobe, Cytocell and Chromoprobe Multiprobe are registered trademarks of Cytocell Ltd. The Chromoprobe principle is covered by international patents WO9314223, EP0623177. The design of the Multiprobe is a registered design and is also covered by a Design Patent No. 420,745. Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

Cytocell Ltd.
4 Technopark
Newmarket Road
Cambridge, CB5 8PB, UK.
T: +44(0)1223 294048
F: +44(0)1223 294986
E: probes@cytoCELL.com
W: www.cytoCELL.com

004/2008-03-20

PI027/CE