Instructions For Use

REF: LPS 003

EGFR

FOR PROFESSIONAL USE ONLY

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
Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) has proven to be a vital adjunct to classical cytogenetics since 1986. More recently, because it can show specific changes in individual cells, it has become useful in the assessment of a vast number of clonal lymphoproliferative disorders as it can light-up DNA sequences in any cell, at any stage of its cell cycle. The process relies on the manipulation of the DNA helix's inherent stability as a double stranded molecule to first separate the strands and then substitute one of them with a synthetic, fluorescently labelled strand of the same sequence. The helix readily rehybridises and thus labels the target sequence with a detectable fluorescent dye. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or Southern blotting can provide information at the gene expression level but, when carried out on tissue sections (either cryostat or paraffin embedded), FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.

Probe Specification

EGFR: Red
7 cen: Green
EGFR (also known as HER1) Epidermal Growth Factor Receptor on chromosome band 7p11.2 is a type I tyrosine kinase receptor for members of the epidermal growth factor family. Binding of the EGFR to epidermal growth factor protein leads to cell proliferation. Abnormally elevated EGFR kinase activity can lead to proliferative diseases such as non-small-cell lung carcinoma (NSCLC), which accounts for 80-85% of all lung cancers¹ and less frequently in breast cancer² among others. There are three small molecule EGFR kinase inhibitor drugs in clinical use (gefitinib, erlotinib and lapatinib), and several others are currently undergoing clinical development³. Approximately, 10% of lung cancer patients have rapid and dramatic response to these TKIs^{4,5}. FISH has proven to be useful for determining the status and therefore selection of NSCLC patients for treatment with EGFR TKIs⁶.

The 295 kb EGFR probe, labelled in red, covers the EGFR gene at 7p11.2 and detects the change in copy number at this region. The centromeric probe in green acts as the control for chromosome 7.

In a normal cell, two red and two green signals are observed while in cells with amplification of the EGFR locus, multiple red signals will be observed.

Materials Provided

Probe: 100µl per vial (10 tests)
Amount of EGFR probe: 160-200 ng/test
Amount of 7cen probe (D7Z1): 240-300 ng/test
The probes are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: EGFR Probe with a red fluorophore (specific to the Texas Red spectrum) and the 7 cen probe (D7Z1) Control Probe with a green fluorophore (specific to the FITC spectrum).
Counterstain: 150 µl per vial (15 tests)
The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 µl to 200 µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendations

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores simultaneously.

Sample Preparation

The kit is designed for use on paraffin embedded tissue sections or tissue microarrays (TMA) and should be prepared according to the laboratory or institution guidelines. For FISH, 4 to 6 µm thick formalin-fixed and paraffin-embedded tissue specimens should be used.

Paraffin Embedded Tissue Sections FISH Suggested Protocol

Paraffin embedded tissue sections require pre-treatment before hybridizing with FISH probes i.e. removal of paraffin wax (de-paraffination), rehydration, acid treatment, pre-treatment and protein digestion.

Reagents:

- > Hydrochloric acid (HCl)
- > 2xSSC (saline sodium citrate) wash buffer
- > HistoClear (e.g. Fisher Scientific (Prod. Code: HIS-010-010S))
- > 70%, 85%, 100% Ethanol
- > Lyophilized pepsin (e.g. SIGMA cat. number: P6887, 3200-4.500 U/mg)
- > Paraffin embedded tissue sections (PETS)
- > Sodium thiocyanate NaSCN (e.g. 81.07g/mol SIGMA cat. number: 251410)
- > 0.4xSSC pH 7.0 wash buffer
- > 2xSSC-0.05% Tween-20 pH 7.0 wash buffer

Reagent Preparation

- Acid Pre-Treatment Solution
 - > Prepare 100ml of 1M NaSCN Pre-treatment solution by dissolving 8.1g in 100ml distilled water
 - > Pour Pre-Treatment solution into 100ml ceramic or glass Coplin jar. Place into water bath and set the temperature at 80°C. Ensure that the temperature of the Pre-Treatment solution is equilibrated to 80°C before commencing the relevant steps.
- 2xSSC Wash Buffer
 - > 2xSSC Buffer should be allowed to reach room temperature before use
- Pepsin
 - > The 100ml Coplin jar of 0.01M HCl should be preheated in 37°C waterbath. Add 0.2g pepsin should be added and mixed **just before** immersing slides in the solution

De-paraffinisation sections (Fume Cabinet)

- > Mark probe target area with a diamond marker pen.
- > Pour 100 ml HistoClear into three ceramic or glass Coplin jars, labeled 616, 626 and 636. Leave the slides in jar 616 for 10 min. Repeat with jars 626 and 636.
- > Remove slides from jar 636 and remove excess solution by gently tapping the edges of each slide on a paper towel.
- > Pour the used HistoClear into waste bottles marked HistoClear 616, 626, 636 when finished. Change the HistoClear after 200 slides have been processed or if cloudy in appearance.

Re-hydration

- > Place slides in a plastic jar containing 50ml 100% ethanol (room temperature) for 2 min. pour off the ethanol into a suitable waste container.
- > Replace with fresh 100% ethanol. After 2 min pour off the ethanol into a suitable waste container.
- > Place slides in a plastic jar containing 50ml of 70% ethanol (room temperature) for 2 min. pour off the ethanol into a suitable waste container.
- > Replace with fresh 70% ethanol. After 2 min pour off the ethanol into a suitable waste container.
- > Air dry slides.

Acid treatment

- > Place slides in 100ml 0.1 M HCl for 20min.
- > Place slides in distilled water for 3 min.
- > Place slides in 2xSSC buffer for 3 min.

Pre-Treatment (Fume Cabinet)

- > Place slides in NaSCN Pre-Treatment solution pre-heated to 80°C in a waterbath and incubate for 30 min.
- > Remove slides from the jar and remove excess solution by gently tapping the edges of the each slide on a paper towel.
- > Place slides in distilled water for 1 min at room temperature (in a plastic jar)
- > Remove slides from the jar and place in 2xSSC Buffer for 3 min at room temperature. Pour off the buffer into a suitable waste container.
- > Replace with fresh 2xSSC buffer. After 3 min. pour of the buffer into a suitable waste container.
- > Remove slides from the jar and remove the excess solution by gently tapping the edges of each slide on a paper towel.

Pepsin Digestion

- > Place slides in the pre-heated pepsin solution and incubate for 10-20 min. (Coplin jar).
- > Remove slides from the jar and gently tap the edges of each slide on a paper towel to remove the Pepsin.
- > Air dry slides.
- > Dehydrate tissue sections through the ethanol series (70%, 85% and 100%) in Coplin jars (room temperature) immersing the slides for 2 min. in each.

Hybridisation

- > Add 10µl of probe on section (the amount of probe depends on a size of a section, add more if needed)
- > Apply a 24x24 coverslip and seal with rubber solution glue.
- > Place the slide on a hot plate and denaturate at 75°C for 5 min.
- > Place the slides in the 37°C incubator over night.

Post Hybridisation Washes

- > Remove coverslip and all traces of glue carefully.
- > Wash slide in 0.4xSSC (pH 7.0) at 72°C for 2 min.
- > Drain slide and wash in 2xSSC/Tween (pH 7.0) at RT for 30 sec.
- > Drain the slide and apply 10-15 1ES DAPI antifade.

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Comments:

Hybridation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridization efficiency (e.g. long incubation times of NaSCN or pepsin digestion) tend to destroy cell structure and tissue morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results. The optimal length of protease treatment time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Protease treatment should be decreased for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. These samples need to be handled particularly carefully to avoid over-digestion.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.