## Instructions For Use

REF: LPS 006

### EWSR1

FOR PROFESSIONAL USE ONLY

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Further information available at [www.cytocell.com](http://www.cytocell.com)


Fluorescence *In Situ* Hybridisation (FISH) has proven to be a vital adjunct to classical cytogenetics since 1986. More recently, because it can show specific changes in individual cells, it has become useful in the assessment of a vast number of clonal lymphoproliferative disorders as it can highlight DNA sequences in any cell, at any stage of its cell cycle. The process relies on the manipulation of the DNA helix's inherent stability as a double stranded molecule to first separate the strands and then substitute one of them with a synthetic, fluorescently labelled strand of the same sequence. The helix readily rehybridises and thus labels the target sequence with a detectable fluorescent dye. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or Southern blotting can provide information at the gene expression level but, when carried out on tissue sections (either cryostat or paraffin embedded), FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.

#### Probe Specification

Ewing's sarcoma is the second most frequent primary bone cancer in patients below 20 years of age<sup>1</sup>. It is characterized mostly by translocation between the two genes EWSR1 (Ewing's Sarcoma Region 1) and FLI1 (Friend leukemia virus integration 1) in the t(11;22)(q24;q12) translocation, found in 85% of such tumours<sup>2</sup>. Ten percent of Ewing's sarcoma cases have a translocation involving EWSR1 and ERG (t(21;22)(q22;q12)), and < 1% of Ewing's sarcomas harbour each of the 7p22 (ETV1), 17q12 (E1AF) and 2q36 (FEV) translocations with EWSR1. EWSR1 is also involved in translocations in desmoplastic small round cell tumour, a subset of angiosarcoma fibrous histiocytomas and possibly in myxoid liposarcoma<sup>3,4</sup>.

The EWSR1 breakapart probe can also be used to distinguish rare soft tissue sarcoma: a clear cell sarcoma from malignant myeloma which is hard to distinguish histologically and immunohistologically. EWS-ATF (Activating transcription factor) translocation t(12;22)(q13;q12), which never happens in the malignant myeloma has been identified in 70% to over 90% of clear cell sarcoma<sup>5</sup>.

FISH has been shown to be a more sensitive and reliable method than RT-PCR for the diagnosis of EWS in solid tissues<sup>6</sup>.

The EWSR1 probe mix contains a green (153 kb) probe and a red (141 kb) probe which are positioned on each side of the EWSR1 gene.

In the normal cell, 2 red/green (can appear as yellow, Y) signals are expected (2Y). A translocation will result in 1R, 1G, 1Y: 1 green and 1 red signal from the translocated chromosomes and 1 fused red/green signals (can appear as yellow) from the normal chromosome.

#### Materials Provided

Probe: 100µl per vial (10 tests)

Amount of EWSR1 probe in red : 80-100 ng/test

Amount of EWSR1 probe in green: 240-300 ng/test

The probes are provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: with a red fluorophore (specific to the Texas Red spectrum) and with a green fluorophore (specific to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

#### Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.

- All hazardous materials should be disposed of according to your institutions' guidelines for hazardous waste disposal.

#### Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

#### Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 µl ó 200 µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

#### Fluorescence Microscope Recommendations

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores simultaneously.

#### Sample Preparation

The kit is designed for use on paraffin embedded tissue sections or tissue microarrays (TMA) and should be prepared according to the laboratory or institution guidelines.

For FISH, 4 to 6 µm thick formalin-fixed and paraffin-embedded tissue specimens should be used.

#### Paraffin Embedded Tissue Sections FISH Suggested Protocol

Paraffin embedded tissue sections require pre-treatment before hybridizing with FISH probes i.e. removal of paraffin wax (de-paraffination), rehydration, acid treatment, pre-treatment and protein digestion.

#### Reagents:

- > Hydrochloric acid (HCl)
- > 2xSSC (saline sodium citrate) wash buffer
- > Histoclear (e.g. Fisher Scientific (Prod. Code: HIS-010-010S))
- > 70%, 85%, 100% Ethanol
- > Lyophilized pepsin (e.g. SIGMA cat. number: P6887, 3200-4.500 U/mg)
- > Paraffin embedded tissue sections (PETS)
- > Sodium thiocyanate NaSCN (e.g. 81.07g/mol SIGMA cat. number:251410)
- > 0.4xSSC pH 7.0 wash buffer
- > 2xSSC-0.05% Tween-20 pH 7.0 wash buffer

#### Reagent Preparation

- Acid Pre-Treatment Solution
  - > Prepare 100ml of 1M NaSCN Pre-treatment solution by dissolving 8.1g in 100ml distilled water
  - > Pour Pre-Treatment solution into 100ml ceramic or glass Coplin jar. Place into water bath and set the temperature at 80°C. Ensure that the temperature of the Pre-Treatment solution is equilibrated to 80°C before commencing the relevant steps.
- 2xSSC Wash Buffer
  - > 2xSSC Buffer should be allowed to reach room temperature before use
- Pepsin
  - > The 100ml Coplin jar of 0.01M HCl should be preheated in 37°C waterbath. Add 0.2g pepsin should be added and mixed **just before** immersing slides in the solution

#### De-paraffinisation sections (Fume Cabinet)

- > Mark probe target area with a diamond marker pen.
- > Pour 100 ml Histoclear into three ceramic or glass Coplin jars, labeled 616, 626 and 636. Leave the slides in jar 616 for 10 min. Repeat with jars 626 and 636.
- > Remove slides from jar 636 and remove excess solution by gently tapping the edges of each slide on a paper towel.
- > Pour the used Histoclear into waste bottles marked Histoclear 616, 626, 636 when finished. Change the Histoclear after 200 slides have been processed or if cloudy in appearance.

#### Re-hydration

- > Place slides in a plastic jar containing 50ml 100% ethanol (room temperature) for 2 min. pour off the ethanol into a suitable waste container.
- > Replace with fresh 100% ethanol. After 2 min pour off the ethanol into a suitable waste container.
- > Place slides in a plastic jar containing 50ml of 70% ethanol (room temperature) for 2 min. pour off the ethanol into a suitable waste container.
- > Replace with fresh 70% ethanol. After 2 min pour off the ethanol into a suitable waste container.
- > Air dry slides.

#### Acid treatment

- > Place slides in 100ml 0.1 M HCl for 20min.
- > Place slides in distilled water for 3 min.
- > Place slides in 2xSSC buffer for 3 min.

#### Pre-Treatment (Fume Cabinet)

- > Place slides in NaSCN Pre-Treatment solution pre-heated to 80°C in a waterbath and incubate for 30 min.
- > Remove slides from the jar and remove excess solution by gently tapping the edges of the each slide on a paper towel.
- > Place slides in distilled water for 1 min at room temperature (in a plastic jar)
- > Remove slides from the jar and place in 2xSSC Buffer for 3 min at room temperature. Pour off the buffer into a suitable waste container.
- > Replace with fresh 2xSSC buffer. After 3 min. pour of the buffer into a suitable waste container.
- > Remove slides from the jar and remove the excess solution by gently by gently tapping the edges of each slide on a paper towel.

#### Pepsin Digestion

- > Place slides in the pre-heated pepsin solution and incubate for 10-20 min. (Coplin jar).
- > Remove slides from the jar and gently tap the edges of each slide on a paper towel to remove the Pepsin.
- > Air dry slides.
- > Dehydrate tissue sections through the ethanol series (70%, 85% and 100%) in Coplin jars (room temperature) immersing the slides for 2 min. in each.

#### Hybridisation

- > Add 10µl of probe on section (the amount of probe depends on a size of a section, add more if needed)
- > Apply a 24x24 coverslip and seal with rubber solution glue.
- > Place the slide on a hot plate and denaturate at 75°C for 5 min.
- > Place the slides in the 37°C incubator overnight.

#### Post Hybridisation Washes

- > Remove coverslip and all traces of glue carefully.
- > Wash slide in 0.4xSSC (pH 7.0) at 72°C for 2 min.

- Drain slide and wash in 2xSSC/Tween (pH 7.0) at RT for 30 sec.
- Drain the slide and apply 10-15 1 ES DAPI antifade.
- Cover with the coverslip.

#### Comments:

Hybridisation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridization efficiency (e.g. long incubation times of NaSCN or pepsin digestion) tend to destroy cell structure and tissue morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results. The optimal length of protease treatment time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Protease treatment should be decreased for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. These samples need to be handled particularly carefully to avoid over-digestion.

#### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

#### Procedural Recommendations

1. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
2. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

#### Customer Support

Please contact the Cytocell Sales and Marketing Department.

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#### Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
**Cytocell Ltd.**