## Instructions For Use

REF: LPU 013

## Cri-Du-Chat and SOTOS Probe Combination

FOR PROFESSIONAL USE ONLY

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Further information available at [www.cytocell.com](http://www.cytocell.com)

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.


### Probe Specification

Cri-Du-Chat : CTNND2 (5p15.2) Red  
FLJ25076 (5p15.31) Green  
SOTOS: NSD1 (5q35) Green

### Introduction

Cri-Du-Chat syndrome consists of multiple congenital anomalies, mental retardation, microcephaly, abnormal face, and a mewing cry in infants. It is associated with deletions, which vary in size, of part of the short arm of chromosome 5<sup>1</sup>. The estimated prevalence varies between 1 in 20,000 to 1 in 50,000 births<sup>2</sup>. A critical chromosomal region involved in the high-pitched cry is mapped to proximal 5p15.3, while the region involved in the remaining features of the syndrome mapped to 5p15.2<sup>3,4,5</sup>.

SOTOS syndrome is a neurological disorder characterised by a distinctive facial appearance,

The FLJ25076 probe is labelled in green and is approximately 194 kb, contains D5S1637E and D5S2678 markers, and covers the entire gene and flanking DNA. The CTNND2 probe is labelled in red, contains the D5S2883 marker and is approximately 159 kb. The SOTOS probe, covering NSD1, is labelled in green and is approximately 167 kb. In the normal cell, there should be fusion of the red and green signals (2Y) and two green signals (2G), whilst a deletion of FLJ25076 probe target results in 1R, 1Y, 2G, a deletion of CTNND2 target results in 1Y, 3G and deletion of NSD1 target results in 2Y, 1G.

### Materials Provided

Probe: 50 µl per vial (5 tests) or 100 µl per vial (10 tests)  
Amount of CTNND2 probe : 86.6 ng/test  
Amount of FLJ25076 probe : 7.3 ng/test  
Amount of NSD1 probe : 86.6 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: CTNND2 Probe with a red fluorophore (specificity to the Texas Red spectrum), the FLJ23076 and NSD1 with a green fluorophore (specificity to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

### Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

### Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

### Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl ó 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

### Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

### Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines

Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

### FISH Protocol

#### Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

#### Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

#### Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

#### Hybridisation:

8. **Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

#### Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

#### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

#### Procedural Recommendations

1. Baking or aging of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

#### Customer Support

Please contact the Cytocell Sales and Marketing Department.

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#### Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
**CytoCELL Ltd.**