Instructions For Use

REF: LPU 014

DiGeorge/TBX1 Probe (22q11) with 22q13.3 Deletion Syndrome Probe Combination

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Specification

TBX1 (22q11.2) Region Probe Red
22qter Subtelomere Specific Probe (N85A3) Probe Green

Introduction


DiGeorge syndrome (DGS)¹ and a variety of congenital malformation syndromes, including Velocardiofacial (VCF)² and Conotruncal Anomaly Face syndromes³ share the phenotypic features covered by the acronym CATCH22 (cardiac defects; abnormal facies; thymic hypoplasia; cleft palate; hypocalcaemia) and deletion of chromosome 22 at 22q11.2^{4,5,6}. In addition around 17% of nonsyndromic patients with isolated conotruncal defects have been shown to have a 22q11.2 microdeletion⁷. The incidence of these anomalies is estimated to be 1:4000 live births⁸ and therefore the deletion 22q11.2 represents one of the most common genetic defects.

A region of 2 Mb referred to as the DiGeorge Critical Region (DGCR) is most commonly deleted in up to 90% of patients^{9,10}. Within the DGCR a minimal critical region of 480-575 kb has been described^{11,12} containing several genes including the T-box transcription factor -TBX1, a key candidate for the aortic arch malformations and anomalies related to pharyngeal apparatus and ear development seen in DGS^{13,14,15}. In addition, TBX1 may play a role in behaviour¹⁶. About 35% patients with DiGeorge develop psychiatric disorders, mainly schizophrenia and bipolar disorder.

22q13.3 deletion syndrome

The 22q13.3 deletion syndrome presents a recognisable phenotype characterised by hypotonia, delay or absence of expressive speech, moderate to profound mental retardation, normal to accelerated growth and mild dysmorphic features¹⁷. Some deletions of the terminal region of chromosome 22q are cytogenetically visible. However, a few cases of cryptic deletions have been reported^{17,18} suggesting that the actual incidence of 22q telomere deletion may be higher than previously thought.

Several observations of patients with 22q13.3 deletion showed that the ProSAP2/SHANK3²⁴ gene, coding for a structural protein of the postsynaptic density of excitation synapses and expressed in the cortex and cerebellum of the brain¹⁹, was disrupted^{19,20,21} or deleted²², making it a good candidate gene for this syndrome. The deletion varies widely in size, from 130kb to 9Mb^{22,23,24}. Recent findings therefore recommend the use of 22q subtelomeric probes distal to the ARSA gene for examining all 22q13.3 deletions^{24,25}.

The TBX1 Region Probe labelled in red, is approximately 213 kb, contains the D22S1627 marker. The 22qter subtelomere specific probe (clone N85A3), labelled in green, is sited in the ProSAP2/SHANK3 gene, allows identification of the most distal 22q13.3 deletions. Both probes allow identification of chromosome 22, acting as control probes. In the normal cell, there should be two red and two green signals (2R, 2G) whilst a deletion of the DGCR probe target results in one red signal and two green controls (1R, 2G) and a deletion of the 22q subtelomeric probe results in two red controls and one green signal (2R, 1G).

Materials Provided

Probe: 50µl per vial (5 tests) or 100 µl per vial (10 tests)

Amount of TBX1 probe : 7.3 ng/test

Amount of 22q probe: 108.2 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: TBX1 Probe with a red fluorophore (specificity to the Texas Red spectrum) and the 22q Control Probe with a green fluorophore (specificity to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen ; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl to 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines

Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

Hybridisation:

8. **Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking or aging of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.