Instructions For Use

REF: LPU 015

DiGeorge II Probe (10p14) with 10 Centromere Control

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Further information available at www.cytocell.com


Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Specification

CUGBP2 (10p14) Region Probe Red
10c Probe Green

Introduction

DiGeorge syndrome¹ and a variety of congenital malformation syndromes, including Velocardiofacial (VCF)² share the phenotypic features covered by the acronym CATCH22 (cardiac defects; abnormal facies; thymic hypoplasia; cleft palate; hypocalcaemia) and deletion of chromosome 22 at 22q11^{2,3,4,5}. However, 10p13.14 deletions have been shown in some patients with DGS^{6,7,8}. The deletion of the DGS2 locus on 10p may be 50 times less frequent than the deletion of the DGS1 locus on 22q and is estimated to occur in 1 in 200,000 live births⁹. A gene called BRUNOL3 has been identified within the 300 kb minimal region of DGS2 and postulated to be involved in DGS2 deletion¹⁰. It is a candidate gene for the heart defect and thymus hypoplasia/aplasia associated with partial monosomy 10p¹⁰. This gene may be involved in atrial septal defects (ASD), which is common cardiac anomaly associated with DGS2¹¹.

The CUGBP2 (BRUNOL3) Region Probe labelled in red, is approximately 144 kb, contains the markers D10S2196 and D10S2420 and covers much of the BRUNOL3 gene. The accompanying 10 centromere Probe (D10Z1), labelled in green, acts as a control probe. In the normal cell, there should be two red and two green signals (2R, 2G) whilst a deletion of a probe target results in one red signal and two green control signals (1R, 2G).

Materials Provided

Probe: 50µl per vial (5 tests) or 100 µl per vial (10 tests)

Amount of CUGBP2 probe : 7.3 ng/test

Amount of 10c probe: 8.7 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: CUGBP2 probe with a red fluorophore (specificity to the Texas Red spectrum) and the 10c Control Probe with a green fluorophore (specificity to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen ; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl to 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

Hybridisation:

8. Hybridise slide overnight in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking or aging of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
- Berend, S.A. *et al* (2000) *Am J Med Genet* 91(4) : 313-317
- Yatsenko, S.A. *et al* (2004) *Clin Genet* 66:128-136
CytoCELL Ltd.