Instructions For Use

REF: LPU 016

Kallmann (KAL1)/ Steroid Sulphatase Deficiency (STS) Probe Combination

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.


Probe Specification

KAL1 (Xp22.3) Red
STS (Xp22.3) Green
Xc: DXZ1 (Xp11.1-q11.1) control probe Green

Introduction

Kallmann syndrome (KS) is a developmental disease characterised by olfactory deficiency and hypogonadotropic hypogandism (HH), which is responsible for the absence of spontaneous puberty¹. It is a heterogeneous developmental genetic disorder affecting about 1 in 8,000 males and 1 in 40,000 females². Reports indicate three modes of inheritance: X-linked, autosomal dominant and autosomal recessive¹³⁴. It has been shown that mutations in KAL1 on Xp22.3 result in the X-linked form⁵. KAL1 consists of 14 exons and extends over 210 kb⁶. KAL1 gene abnormalities reported in patients with KS include missense and nonsense mutations, splice site mutations, intragenic deletions and submicroscopic chromosomal deletions involving the entire KAL1 gene⁷.

Steroid Sulphatase Deficiency (STS) (also known as X-linked Ichthyosis)⁸ is the second most common type of ichthyosis and one of the most frequent human enzyme deficiency disorders. Deficiency of the STS enzyme⁹ is known to be responsible for dark, adhesive and regular scaling of the skin. The STS gene is mapped at the distal part of the short arm of the X chromosome, a region which escapes X-chromosome inactivation and has the highest ratio of chromosomal deletions among all genetic disorders¹⁰. Complete deletions have been found in more than 90% of patients¹¹. The deletions can extend to involve neighbouring genes causing contiguous gene defect. Therefore, STS may be associated with KS¹².

The KAL1 Region probe labelled in red, is approximately 336 kb, it contains markers DXS278, DXS7470, DXS7714, DXS7053 and covers the entire KAL1 gene and flanking DNA. The accompanying STS Region probe, labelled in green, consists of two clones: a 173 kb clone containing marker DXS7731 and a 149 kb clone containing markers DXS6767, DXS7965, DXS7500. They cover the entire STS gene and flanking DNA. The X centromere probe: DXZ1 (Xp11.1-q11.1) probe, labelled in green, acts as a control. In the normal XY cell, there should be one red and two green signals (1R, 2G) or fusion of the red and green signals (yellow) and one green control signal (1Y, 1G), whilst a deletion of a probe target results in no red/yellow signal and one green control signal (0R, 1G) (KAL1/STS deletion) or (1R, 1G)(STS deletion) or (0R, 2G) (KAL1 deletion). In the normal XX cell, there should be two red and four green (2R, 4G) or two yellow and two green (2Y, 2G) signals, whilst a deletion of a probe target results in 1Y, 2G (KAL1/STS deletion) or 1Y, 1R, 2G (STS deletion) or 1Y, 3G (KAL1 deletion).

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)
Amount of KAL1 probe: 12.2 ng/test
Amount of STS probe: 144.2 ng/test
Amount of Xc probe: 8.7 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: KAL1 Probe with a red fluorophore (specificity to the Texas Red spectrum) and the STS and Xc Control Probe with a green fluorophore (specificity to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl to 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

Hybridisation:

8. **Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking and aging of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Comment

The KAL1 region probe may show faint cross-hybridisation to the Y chromosome, because of sequence homologies in this region.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.