Instructions For Use

REF: LPU 019

Smith-Magenis (RAI1) and Miller-Dieker (LIS1) Probe Combination

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.


Probe Probe Specification
RAI1 (SMCR 17p11.2) Green
LIS1 (MDS/ILS 17p13.3) Red

Introduction

Smith-Magenis syndrome (SMS) is a multiple congenital anomaly syndrome characterised by mental retardation, neuro-behavioral abnormalities, sleep disturbances, short stature, minor craniofacial and skeletal anomalies, congenital heart defects and renal anomalies^{1,2}. It is one of the most frequently observed human microdeletion syndromes and associated with an interstitial deletion of the chromosome band 17p11.2^{1,2}. Molecular studies in SMS patients suggest a common deletion region spanning approximately 700 kb³. The proximal boundary is within a region of overlap between the FLI1 and LLLGL1 genes and the distal boundary within the PEMT gene³. Deletions or mutations in RAI1 (retinoic acid induced 1) gene, which lies within the 17p11.2 locus were associated with the syndrome^{3,4,5,6}. RAI1 was shown to be the primary gene responsible for most features of SMS^{7,8}.

Miller-Dieker syndrome (MDS) is a multiple malformation characterised by classical lissencephaly, a characteristic facial appearance and sometimes other birth defects⁹. It is associated with visible or submicroscopic rearrangements within chromosome band 17p13.3 in almost all cases¹⁰. Isolated lissencephaly sequence (ILS) consists of classical lissencephaly with no other major anomalies¹¹. Submicroscopic deletions of chromosome 17p13.3 have been detected in almost 40% of these patients¹⁰.

MDS is considered a contiguous gene deletion syndrome where deletion of physically contiguous genes leads to the complex phenotypic abnormalities seen in MDS. LIS1 is located at 17p13.3 and recognised as the causative gene responsible for the lissencephaly phenotype in both MDS and ILS^{12,13}. A deletion of MDS patients always involves LIS1 together with telomeric loci to in excess of 250 kb¹².

The SMCR probe is 160 kb and targets the RAI1 gene including the D17S258 marker. The MDS/ILS probe is 90 kb and targets the entire LIS1 gene of 80 kb. This probe product is dual labelled and intended to identify deletions of SMS and deletions or parental cryptic translocations of MDS/ILS patients. The two unique sequences provide each other with a control and allow identification of chromosome 17. In the normal cell, there should be two red and two green signals (2R, 2G), whilst a deletion of a probe target results in 1R, 2G (MDS) or 2R, 1G (SMS) or 1R1G.

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of RAI1 (SMCR) probe: 86.6 ng/test

Amount of LIS1 (MDS/ILS) probe: 7.3 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: RAI1 (SMCR) with a green fluorophore

(specificity to the FITC spectrum) and LIS1 (MDS/ILS) with red fluorophore (specificity to the Texas Red spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl to 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines. Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

Hybridisation:

8. **Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking or aging of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.