Instructions For Use

REF: LPU 023

Rubinstein-Taybi

FOR PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com


Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Specification

CREBBP (16p13.3) Red
16 cen (D16Z2) Green

Introduction

Rubinstein-Taybi syndrome described in 1963 is a well-defined condition characterized by short stature, mental retardation, facial abnormalities, and broad thumbs and first toes¹. This condition is an autosomal-dominant disease with a frequency of one individual in 100,000 newborns. Chromosome band 16p13.3 contains the Rubinstein-Taybi syndrome locus, which includes a gene encoding a binding protein for cyclic adenosine monophosphate response element binding protein (CBP) (CREBBP). This gene spans approximately 150 kb with 31 exons² and dosage anomalies are the common cause of the disease³. Disruptions of CREBBP, either by gross chromosomal rearrangement or point mutation, are responsible for the phenotype of Rubinstein-Taybi syndrome. Microdeletions including the CREBBP gene account for about 10-15% of RTS cases¹.

EnsEMBL 54

The CREBBP region probe labelled in red, is approximately 300 kb and covers the entire CREBBP gene and flanking DNA. The accompanying D16Z2 centromere specific probe in green, acts as a control probe. In the normal cell there should be two red and two green signals (2R, 2G), whilst a deletion of a probe target results in one red signal and two green control signals (1R, 2G).

Materials Provided

Probe: 50µl per vial (5 tests) or 100 µl per vial (10 tests)

Amount of CREBBP probe : 40-50 ng/test

Amount of 16 cen probe: 120-150 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: CREBBP Probe with a red fluorophore (specificity to the Texas Red spectrum) and the 16 cen Specific Control Probe with a green fluorophore (specificity to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.

- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C)
- Variable volume micropipettes range 1 µl to 200 µl
- Water bath with accurate temperature control at 72°C
- Microcentrifuge tubes (0.5 ml)
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- Plastic or glass coplin jars
- Forceps
- Fluorescence grade microscope lens immersion oil
- Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

Sample Preparation

The probe is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines.

Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

- Spot cell sample onto cleaned microscope slide.
- Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation

- Remove probe from -20°C freezer and allow to warm to RT.
- Ensure probe solution is uniform by repeated pipette mixing.
- Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

Denaturation

- Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

Hybridisation:

- Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes

- Remove coverslip and all traces of glue carefully.
- Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
- Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
- Drain the slide and apply 10 µl of DAPI antifade.
- Cover with a coverslip and allow colour to develop in the dark for 10 mins.
- View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

- Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
- The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.