## Instructions For Use

REF: LPU 024

## Saethre-Chotzen / Williams-Beuren Combination

FOR PROFESSIONAL USE ONLY

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Further information available at [www.cytocell.com](http://www.cytocell.com)

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.


### Probe Specification

TWIST1 (7p21.1) Red  
WBS CR (7q11.23) Green

### Introduction

One of the most frequent syndromes in congenital malformation: craniosynostosis (premature closure of one or more cranial sutures) is Saethre-Chotzen syndrome (acrocephalosyndactyly type III; ACS III). Saethre-Chotzen syndrome is characterized by craniofacial and limb abnormalities (Saethre M 1931, Chotzen F 1932) and the incidence of this rare syndrome is estimated at between 1 in 25-50,000 live births. This autosomal dominant disorder is one of the most difficult to diagnose and the identification of the TWIST1, a basic helix-loop-helix transcription factor on chromosome band 7p21.1, as a causative gene has proved invaluable for the diagnosis<sup>1,2</sup>.

Williams-Beuren Syndrome (WBS) is a developmental disorder caused by a microdeletion within band 7q11.23<sup>3</sup>. Patients display connective tissue problems, typically supravalvular aortic stenosis (SVAS), growth retardation, renal anomalies, transient hypercalcaemia, hyperacusis and mental retardation<sup>4</sup>. Haploinsufficiency of the elastin (ELN) gene has been identified as being responsible for SVAS<sup>5,6</sup> however, of the remaining 15 genes identified within the WBS deletion one of the remaining clinical features have been conclusively attributed to any one of these genes. The most recent work maps the WBS deletion region to be 1.6 Mb<sup>7</sup> and flanked by highly homologous duplications of 320-500 kb within which the common breakpoints cluster. The common WBS deletion results from non-homologous recombination between the GTF2I/NCF1 locus and the GTF2IP1/NCF1P1 locus or rare intrachromosomal exchange between the centromeric GTF2IP2/NCF1 and the telomeric GTF2IP2/NCF1P2 duplications.

The TWIST1 region probe labelled in red, is approximately 181 kb and covers the entire TWIST1 gene and flanking DNA. The accompanying Williams-Beuren region probe, labelled in green is approximately 450 kb and consists of three non-overlapping clones which cover much of the deletion region. In the normal cell there should be two red and two green signals (2R, 2G), whilst a deletion of a probe target results in either 1R 2G when TWIST1 is deleted, or 2R 1G when Williams-Beuren region is deleted.

### Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of TWIST1 probe : 95-120 ng/test

Amount of Williams probe: 105-130 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: TWIST1 Probe with a red fluorophore (specificity to the Texas Red spectrum) and the Williams Probe with a green fluorophore (specificity to the FITC spectrum).

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

### Warnings and Precautions

1. Wear gloves when handling DNA probes and DAPI counterstain.
2. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
3. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
4. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

### Storage and Handling

The probe should be stored at 620°C until the expiry date indicated on the label. The probe and counterstain vials must be stored in the dark.

### Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl ó 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

### Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously.

### Sample Preparation

The probe is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines.

Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

### FISH Protocol

#### Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

#### Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

#### Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

#### Hybridisation:

8. **Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

#### Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

#### Stability of Finished Slides

FISH slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

#### Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

#### Customer Support

Please contact the Cytocell Sales and Marketing Department.

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#### Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
**CytoCELL Ltd.**