Instructions For Use

REF: LPU 025

SHOX

FOR PROFESSIONAL USE ONLY

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
Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation and denaturation is available for annealing to a similarly denatured, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of rapid formamide-free stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Specification
SHOX (Xp22.33/Yp11.2) Red
Yqh (DYZ1) Green
X cen (DXZ1) Blue

Introduction
The SHOX (Short stature Homeobox-containing gene) is positioned at Xp22 and Yp11.3 in the pseudoautosomal region (PAR1) of both chromosomes¹. It encodes a transcription factor of 293 and 225 amino acids (SHOXa and SHOXb respectively, whose proteins differ in the C-terminal region). SHOX is a cell-specific homeodomain protein involved in cell cycle and growth regulation and activates transcription in osteogenic cells². SHOX haploinsufficiency is involved in the etiology of idiopathic short stature and the short stature in Turner syndrome while homozygous loss of the SHOX gene has been correlated with the Langer type mesomorphic dysplasia. Subsequently heterozygous SHOX mutations were also shown to cause Leri-Wassilow dyschondrosteosis. The incidence of SHOX deficiency is between 1 in 2000 to 1 in 5000 in the general population and 1 in 40 to 1 in 150 in short people^{3,4}.

EnrEMBL 54

The SHOX region probe labelled in red, is approximately 166 kb and covers the entire SHOX gene and flanking DNA on chromosomes Y and X. The accompanying DYZ1 Y heterochromatin specific probe in green on the Y chromosome and DXZ1 centomere specific probe in blue on the X chromosome, act as the chromosome identification control probes. In the normal male cell there should be two red and one green and one blue signals (2R, 1G, 1B), whilst a deletion of a probe target results in one red signal, one green and one blue control signals (1R, 1G, 1B). In the normal female cell there should be two red and two blue signals (2R, 2B), whilst a deletion of a probe target results in one red signal, two blue control signals (1R, 2B).

Materials Provided

Probe: 50µl per vial (5 tests) or 100 µl per vial (10 tests)

Amount of SHOX probe : 40-50 ng/test

Amount of Yqh probe : 60-75 ng/test

Amount of Xc probe: 120-150 ng/test

The probe is provided premixed and ready to use in hybridisation solution (Formamide; Dextran Sulphate; SSC). The probe DNA is directly labelled: SHOX Probe with a red fluorophore (specificity

to the Texas Red spectrum), the Yqh Probe with a green fluorophore (specificity to the FITC spectrum) and the Xc Probe with a blue fluorophore (specificity to the Aqua spectrum)

Counterstain: 150 µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius kit should be stored at 620°C until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

- a) Hotplate (with a solid plate and accurate temperature control up to 80°C)
- b) Variable volume micropipettes range 1 µl to 200 µl
- c) Water bath with accurate temperature control at 72°C
- d) Microcentrifuge tubes (0.5 ml)
- e) Fluorescence microscope (Please see Fluorescence Microscope Recommendation section)
- f) Plastic or glass coplin jars
- g) Forceps
- h) Fluorescence grade microscope lens immersion oil
- i) Bench top centrifuge

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100 watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all three fluorophores simultaneously. Probe's blue fluorophore has specificity to the Aqua spectrum (Aqua filter is required).

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative and should be prepared according to the laboratory or institution guidelines.

Prepare blood metaphase spreads on microscope slides according to standard cytogenetic procedures.

FISH Protocol

Slide preparation

1. Spot cell sample onto cleaned microscope slide.
2. Immerse slide in 2 x SSC, pH 7.0 for 2 mins.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 mins.

Pre-Denaturation

4. Remove probe from -20°C freezer and allow to warm to RT.
5. Ensure probe solution is uniform by repeated pipette mixing.
6. Remove 10 µl of probe and place on cell sample slide, cover with a 24 x 24mm glass coverslip and seal with rubber solution glue.

Denaturation

7. Place slide on to a 75°C (+/-1°C) hotplate and denature for 2 minutes.

Hybridisation:

8. **Hybridise slide overnight** in a humid, lightproof container at 37°C (+/-1°C).

Post-Hybridisation Washes

9. Remove coverslip and all traces of glue carefully.
10. Wash slide in 0.4 x SSC at 72°C (+/-1°C) (pH 7.0) for 2 mins.
11. Drain slide and wash in 2 x SSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds.
12. Drain the slide and apply 10 µl of DAPI antifade.
13. Cover with a coverslip and allow colour to develop in the dark for 10 mins.
14. View with fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
3. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

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Patents and Trademarks
Aquarius and CytoCELL are registered trademarks of CytoCELL Ltd.
CytoCELL Ltd.